
DUCAGEL Z1

CE

UGIN
dentaire

14 -19 rue Henri-Dunant 38180 SEYSSINS – France
Tél. National : 04.76.84.45.45 - Fax : 04.76.84.46.46
Tél. International : +33.4.76.84.45.43 - Fax : +33.4.76.84.46.46
Email : info@ugin-dentaire.fr - [http : www.ugin-dentaire.fr](http://www.ugin-dentaire.fr)

DUCAGEL Z1

MODE D'EMPLOI	1
INSTRUCTIONS FOR USE	7
GEBRAUCHSANWEISUNG	13
ISTRUZIONI D'USO	19
MODO DE EMPLEO	25
Правила технической эксплуатации	31

MODE D'EMPLOI

SOMMAIRE

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INSTRUCTIONS FOR USE

CONTENTS

INSTALLATION/SPECIFICATIONS:	8
CONTROL PANEL:	9
STARTING AND PROGRAMMING:	10
OPERATION:	11
MAINTENANCE AND RECOMMENDATIONS:	12

INSTALLATION

1. Check that the machine is in perfect condition immediately after unpacking and notify the carrier and supplier of any damage that may have occurred. Indicate any reserves on the consignment note and have the note signed by the carrier; confirm the reserves to the carrier in writing by recorded delivery letter.
2. Position the machine in a well ventilated area with a gap of at least 10 cm between it and any wall.
3. Plug the power lead into the back of the machine.
4. Make sure that the mains voltage matches the voltage on the machine rating plate and then plug the machine into a standard 10/16A socket that includes an earth terminal.

IMPORTANT

The "DUCAGEL Z1" gelatin melter is designed for use in dental prostheses laboratories. Improper use will render the manufacturer's warranty null and void.

SPECIFICATIONS

- Supply voltage: 230 v 50/60Hz
- Power: 850W
- Max. temperature in the tank: 99°C

Dimensions

Height (closed lid):	510mm
Height (open lid):	820mm
Width:	325mm
Height:	425mm

Tank dimensions

Inside diameter:	240mm
Height:	220mm
Capacity:	6 kg of gelatin

Weight: 24.5 kg

-The **Ducagel Z1** has a mixing speed of 30 rpm and includes electronic torque control to prevent damage to the drive unit.


-10 Programs:

N^R. 0 Continuous mixing action during the low-temperature holding level.

N^R. 1 to 9 Intermittent mixing action.

- Delayed start time up to 99 h

CONTROL PANEL:


1. Programme number display (0 to 9).
2. Temperature/Time display.
3. Decrement key.
4. Increment key.
5. Programme start.
6. Stabilisation at low temperature (b).
7. Stop program or current phase.
8. Access to programming function.
9. Drive unit "On" indicator.
10. Fan "On" indicator.
11. Heater "On" indicator.

STARTING

- Fill the tank with gelatin (max. 6 kg) and turn the On/Off switch at the back of the machine to On.
- The temperature in the tank and the program number will appear after a few seconds.
- The programs can now be configured.


PROGRAMMING

- Use the  and  keys to select the required program.
- Press the  key once to set the delayed cycle start (the letter "d". will appear on the display 1 and 0.0 on the display 2).
- The delayed start time setting can be adjusted up to 99h00'.
To be efficient, it's better if its minimum value is not less than 60'.
Use the  and  keys to set the delayed start time.


PLEASE NOTE:

- . If the programmed setting is 0, the cycle will start as soon as the  key is pressed.
- . For a delayed start time, press the  key and the countdown will start before the beginning of the required cycle.

- Press the  key a second time to set the delayed cycle start (the letter "b". will appear on the display 1).
- The low stabilisation temperature setting can be adjusted from 30°C to 75°C (preset to 45°C). Use the  and  keys to alter the setting.
- Press the  key a third time to set the high temperature ("h" will appear at the left of the window).
- The high temperature setting can be adjusted from 75°C to 99°C (preset to 94°C).
Use the  and  keys to set the required temperature.

- Press  the key a fourth time to set the stabilisation time of the top level ("t." will appear at the left of the window).
The stabilisation time can be set to between 10' and 30' (preset to 10') using the  and  keys.
- The programming procedure is now complete. Press the  key to confirm the programmed settings.
- Repeat as above for each program.

OPERATION

- Select a program and then press the  key to run the cycle.
- The heating element indicator light will come on and the heating process will begin.
- The drive unit indicator light  will come on when the temperature in the tank reaches 35°C, and the blade will start to turn intermittently.
- The temperature will remain constant when it reaches the upper set point (**h**) until the end of the stabilisation time (**t**).
- At the end of the stabilisation time, the fan will cut in and the machine will continue to mix with the heating cycle off.
- The fan will stop when the temperature has dropped to the low set point (**b**) and the drive unit will operate intermittently while the product is being stabilised at the low temperature (except for the program **0**).

Caution! For safety reasons, the cycle will not start if the lid is open. In this case "OPEN" will appear on the display, the motor and the heating process are "Off".

PLEASE NOTE

- A dot at the bottom right of the program number indicates that a stabilisation phase is running.
- During program **N^R. 0** the blade will rotate continuously until the end of the low holding level.

IMPORTANT

If necessary, you can switch directly to the low-temperature stabilisation phase (b) by pressing the  key followed by the  key.

⇒ Press  again and the  key to return to the start of the cycle.

⇒ Pressing  interrupts the current program and shuts off the drive unit, fan and heating cycle causing the temperature to drop naturally.

MAINTENANCE

⇒ The **DUCAGEL Z1** machine requires no special maintenance apart from external cleaning (do not use abrasive products).

RECOMMENDATIONS

⇒ Break the gelatin down into small pieces so that the drive unit can restart more efficiently and to reduce the time required to reach liquefaction temperature.

⇒ After a long period of stabilisation, a small amount of water should be added to compensate for any natural evaporation.

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