# DUCATRON SERIE3

CE

**UGIN'**  
dentaire

14 -19 rue Henri-Dunant 38180 SEYSSINS – France

Tél. National : 04.76.84.45.45 - Fax : 04.76.84.46.46

Tél. International : +33.4.76.84.45.43 - Fax : +33.4.76.84.46.46

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# **INSTRUCTIONS FOR USE**

<b>INTRODUCTION</b>	<b>9</b>
<b>TECHNICAL CHARACTERISTICS</b>	<b>9</b>
<b>INSTALLATION AND GETTING STARTED</b>	<b>10</b>
<b>DESCRIPTION OF CONTROL STATION</b>	<b>11</b>
<b>TRIAL RUN(unloaded)</b>	<b>12</b>
<b>INDUCTION GENERATOR</b>	<b>12</b>
<b>ARM BALANCING</b>	<b>13</b>
<b>CHOOSING THE CRUCIBLE</b>	<b>13</b>
<b>MELTING AND INJECTION</b>	<b>14</b>
<b>ARGON GAS</b>	<b>14</b>

## INTRODUCTION

We would like to congratulate you on becoming the new owner of a “Ducatron serie 3”. This machine has been carefully designed, manufactured and tested to give you complete satisfaction.

You will be making perfect castings and risk free, whatever alloy you use (except titanium).

To make the best use of your Ducatron, please read the operating instructions and advice contained in this manual attentively.

## TECHNICAL CHARACTERISTICS

- Power supply : 230 V single -phase
- Frequency : 50/60 Hz
- Power consumption : 2500 W
- Generator operating frequency : 135 kHz
- Maximum temperature : 1500 / 1600 °C
- Arm rotation speed: 500 rpm
- Crucible capacity : 60 g precious or non precious alloy
- Built – in variable power control
- Argon gas: protects the alloy from oxidation.
- **cooling :**
  - c) Electronic generator :  
One 120 x 120 screw fan, output 25 L/s.
  - d) Inductor :  
Close –circuit water pump with 18 litres tank.

### SIZE AND WEIGHT

Width	690 mm
Depth	585 mm
Height w/cover down	950 mm
Height w/cover up	982 mm
Net weight	110 kg


### ACCESSORIES

2 Grey "silicon carbide" crucibles	6 Cradles (6 sizes)
3 White "silicon dioxide" crucibles	1 Set of casting rings
1 Vitrified carbon sheath	1 Conveyor fork
1 Quartz rod	1 Set of fuses
1 Plug	1 Set of adjustable feet
1 Argon joint	

# INSTALLATION AND GETTING STARTED

- 1) Unpack the machine and check that it is in perfect condition; notify the shipper and supplier of any damage. Note any reserves on the shipping order, which should be countersigned by the deliveryman and confirm them to the shipper by registered mail with acknowledgement of receipt.
- 2) Set the machine on a flat surface and make sure it is levelled.
- 3) Unscrew the 4 screws on the panel. Remove the panel
- 4) Take out the water pump located above the grey plastic tank after disconnecting it from the electric plug. Fill the 18 litres tank with water. Put the water pump back in place and close the front panel. Cover the screw holes with the plastic plugs.
- 5) Raise the cover and remove the arm fastener.
- 6) After checking that the mains voltage is the same as that marked on the black plate, plug the cord into a standard 10/16 A grounded wall outlet.

## DESCRIPTION OF CONTROL STATION

- 1 – Power switch **on/off**
- 2 – Arm speed selectors.
- 3 – ARM: **AXIS** indicator.
- 4 – Control button to raise **the induction coil**.  
Function only when arm is centred.
- 5 – Power control dial for induction generator.
- 6 – **Melt** indicator.
- 7 – **Melt** indicator.
- 8 – **Injection** button.  
Function only when cover is closed.
- 9 – **Injection** indicator.
- 10 – **Stop** button.

## TRIAL RUN (unloaded)

Before you melt, we recommend that you test the machine in operation .The following steps should be followed.

- 1) Put the stainless steel ring (Ø 30) into the small cradle.
- 2) Use the fork to position the cradle + ring assembly on the arm.
- 3) Place the silica crucible with the vitrified carbon sheath on the crucible tray.
- 4) Turn the switch to « **ON** » (**Marche**)
- 5) Centre the arm (a signal sounds when it is correctly positioned)
- 6) Press the **SELF** button
- 7) Press the **FUSION** button
- 8) Set the power variation dial to 8.
- 9) Let the crucible temperature rise for 1 minute.
- 10)Close the cover.
- 11)Press the **INJECT** button.
- 12)Let the machine run for 10 to 20 seconds.
- 13)Press the **STOP** button.
- 14)Wait15 seconds (safety margin following rotation).
- 15)Open the cover.

## INDUCTION GENERATOR

Induction heating is based on an entirely different principle from resistance heating. The inductor generates a high-frequency magnetic field which acts directly on the metals to be melted.

The current passing through the metal charge heats and finally melts it.


The generator is set at the right operation frequency to mix the molten metals thoroughly.

Electromagnetic mixing renders a homogeneous alloy and prevents any of the surface overheating characteristic of conventional generators.

The Ducatron Serie 3 generator is equipped with the most advanced semiconductor power gates (Power –MOS).

It has a particularly high efficiency, equivalent to 90% of absorbed power.

It consists of two different parts.

## BALANCING THE ARM

The balancing system consists of a counterweight mounted on the far end of the arm. On the other end, the balance is established with «ring +cradle» which has a constant Weight (1150±70g).

The chart below gives the theoretical weights of the rings (after a run through the furnace) and corresponding cradles.

DIMENSIONS in MM		RING WEIGHT in GRS		CRADLE WEIGHT in Gr
Ring Diam	Height	Metal	No metal	
30	55	100	54	1074
40	55	170	96	1022
50	55	230	150	966
65	55	355	255	843
80	55	500	385	700
90	70		623	527

### IMPORTANT

Never allow the arm to rotate without the «cradle +ring» assembly.

## CHOOSING THE CRUCIBLE

The Ducatron Serie 3 is delivered with two types of crucibles and a vitrified carbon sheath.

### A – SILICON CARBIDE CRUCIBLE (grey)

The silicon carbide crucible is suitable for casting base metals and precious metals

It has high heat-resistant properties for use at temperatures up to 1,500°C or 1,600°C.

The silicon carbide prevents molten metals from sticking between the crucible bottom and the injection nozzle.

### B – SILICON DIOXIDE CRUCIBLE (white)

The silicon dioxide crucible is suitable for casting base metals and precious metals.

It is slightly less temperature resistant (1400°C to 1500°C) than the silicon carbide crucible, but less expensive.

### C – VITREOUS CARBON SHEATH

The vitreous carbon sheath fits inside silicon dioxide and silicon carbide crucibles to protect precious and semi-precious alloys during the melting process and to help improve temperature consistency.

It is non-porous and so protects the alloy from contamination by the carbon.

### NOTA :

The vitrified carbon sheath should never be used for non precious alloys.

## MELTING AND INJECTION

- Flip the min switch to «on» (Marche)
- Select the arm start-up speed.
- Put the metal to be melted in the bottom of the crucible.
- Place the crucible on the insulating tray located on the arm.
- Position the “ring +cradle “assembly on the arm.
- Centre the arm and check that a signal confirms proper positioning.
- Raise the induction coil.
- Close the cover.
- Begin the melting operation.

The mean melting time is 1 minute for 25g of metal, with the vitreous carbon sheath, except BERYLLIUM alloy.

When electromagnetic forces begin to melt the metal, the ingots group together and form a ball. This is the ideal time to begin injection.

- Stop centrifuging after 10 to 20 seconds by pressing **STOP**
- Wait until the cover safety latch is released (20s)
- Open the cover.

## ARGON GAS

### A – FUNCTION

Argon gas is injected during the operation to avoid oxidation of the alloy.

We recommend that you contact “AIR LIQUIDE COMPAGNY” for the argon gas. They will supply you with the right cylinder as well as a flow meter.

### Characteristics

- Cylinders are leased for 2 m<sup>3</sup>, 4 m<sup>3</sup> and more.
- The gas used is: Nerthal, an industrial argon.
- Flow meter : class 0-15 bars.

### B – METHOD

Connect the flow meter to the receptacle at the back of the machine.

Open the cylinder.

Set the flow meter to 5 litres /minute.

Argon gas will automatically be injected as soon as melting begins.

When you'll have finished the casting, shut off the cylinder.

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