



Since 1993

ANAND MEDIPRODUCTS PVT LTD (AMPRO)

AMPRO is a group company of the ANAND GROUP. ANAND Group is actively engaged in a vast range of products including Dental equipment and materials. ANAND MEDICAIDS private limited, one of the group company, is the leader in manufacturing and export of SUCTION EQUIPMENT since 1969. The Brand 'ANAND' is popular with all major hospitals in India and is exported to more than 30 countries worldwide. Other group activities can be seen our websites www.anandind.com, www.suction-india.com, www.dentallabindia.com.

AMPRO was started in 1993 with the sole aim of providing WORLD CLASS DENTAL EQUIPMENT to the Dental fraternity of India. In the 17 years of inception AMPRO has set up a number of Dental Colleges and Dental Laboratories all over India.

AMPRO specializes in designing & equipping dental colleges & laboratories on turn-key basis. Our consultants are qualified to plan & design the colleges as per DCI (Dental Council of India) regulations.

AMPRO has developed a wide range of laboratory products & phantom head simulator units for dental colleges.

AMPRO undertakes regular ceramic training courses in India & abroad for updating its clients on the latest developments in dental laboratory techniques.

Contact your nearest Dental Depot or write to us at info@amprodental.com for your requirements of a AMPRO world-class product.

For more information, kindly visit us at www.amprodental.com

CONTENTS

Casting Machines _____	2	Steam Cleaner _____	13
Ducatron Quattro , Induction Casting Machine		Pressojet	
Minicast , Bench Top Casting Machine		Gel Duplicating Machine _____	13
Centrifugal Casting Machine , Lever Type		Ducagel - ZI	
Ceramic Furnaces & Vacuum Pump ____	3-4	Hydraulic Flask Press _____	13
Satellis , Ceramic & Pressing Furnace		Electrolytic Polishing Unit _____	14
Artis , Furnace for firing high and low fusing ceramics		Vibrator _____	14
Ellipse , Ceramic Furnace		Magnifying Lamp _____	14
FV Plus , Ceramic Furnace		Bench Lighting	
Vacuum Pump , Oil-free Piston Pump		Plaster Dispenser _____	14
Burnout Furnaces _____	5	Phantom Head / Teaching Aids _____	15
Programix 50 , Burnout / Pre-heating Furnace		Phantom Head Simulator Unit	
Burnout Furnace		Phantom Head , Table Model	
Vacuum Mixers _____	6	Jaw Set , With Plastic Articulator	
Ugivac , Vacuum Mixer		Jaw Set , With Metal Articulator	
Mix-R , Vacuum Mixer		Dental Furniture: Clinic / Lab _____	16
Vacuum Mixer		Linea 1 , Workstaton For Clinic	
Sand Blasters / Accessories _____	7	Linea 3 , Workstaton For Clinic	
W.A.F.I.S. , Water Air Filtering System		Sole 1 , Workstation For Lab	
Base 3 Sand Blaster , With two tanks with recycling nozzle		Suction Units _____	17
Sand Blaster		EUROVAC , Portable (Single jar)	
Suctions for Sandblaster _____	8	MB-36 , Heavy Duty (Double jar)	
Pro 3 Shake , Suction for Sandblaster		Alumina / Zirconia / Alloy Coping _____	18
Pro 3 , Suction for Sand Blaster		WOL-CERAM ONE® , Electrolayering CNC Centre	
Surveyor / Milling _____	8	WOL-TEMP 1200 M , Sintering Furnace	
Surveyors / Parallelometers , Surveyors "B" Type		Consumables _____	19
ISO Color Top , Surveyor & Milling Unit		Ceramic - Performance Plus Press	
Model Trimmers _____	9-10	Investment - Ugivest 38	
Model Trimmer Dry , Single Disk		Investment - Ugicast 24	
Model Trimmer Wet , Single Disk		Alumina Sand	
Model Trimmer Dry/Wet , Double Disk		Crucible	
Palatial / Arch Trimmer		Wax/Ugifil	
Lathe / Polishers _____	11	Ugirex III	
High Speed Alloy Grinder		Ring liner	
Aspiration Hood , Polishing Lathe with suction		Turnkey Dental College / Lab Set Up ____	20
Dental Lathe		Client List _____	21
Dental Hanging Engine			
Laser Pin / Die Cutting _____	12		
Laser Pin , Pindex Machine Laser Pin-Drilling Unit			
Laser Pin , AMPRO			
Die Cutting Unit , Trunk Separator			



« Ducatron Quattro Induction Casting Machine

MAIN FEATURES:

- Higher starting torque.
- Closed circuit cooling system. The coolant level can be checked annually using the front mounted operating panel.
- Coolant level warning light.
- Cover lock indicator.
- Clear view of the alloy fusion process.
- Touch button control panel.
- Timed end-of-cycle shutdown of the cooling pump.
- Auto standby function when idle for 5 minutes.
- Standard accessories included

TECHNICAL DATA

Power Supply	: 230 V - 50/60 Hz, 2500 W
Generator Power Consumption	: 2000 W
Dimension	: 95 kg
Weight	: 600 x 980 x 620 mm



NEW UPGRADED MODEL



Minicast Bench Top Casting Machine

MAIN FEATURES:

- The alloy is melted by using a blowtorch (supplied as an option).
- The casting arm is driven by an electric motor.
- The rotation time can be programmed.
- The cycle is automatic.
- The protective cover is locked during the centrifuging time.

TECHNICAL DATA

Rotation Speed	: 500 rpm/min
Rotation Time	: 6 to 60 s
Power	: 230 V- 50/60Hz, 400 W
Weight	: 50 kg
Dimension	: 330 x 700 x 530 mm



AMPRO Centrifugal Casting Machine



MAIN FEATURES:

- Releasing lever for safety
- Four stainless steel flask rings
- 2 crucible attachments with crucibles one for small crucible and one for big crucible
- 4 rubber ring base former
- One aluminium ring base former
- 3 stainless steel cradles for holding the ring with one fixed balance weight

TECHNICAL DATA

Height	: 26 cms
Die of Base Plate	: 17 cms
Lever Length	: 10"
Arm	: 39 cms
Power Supply	: 230 V- 50/60 Hz
Weight	: 7.5 kg
Dimension	: 26 x 17 x 39cm



« Satellis Ceramic & Pressing Furnace

MAIN FEATURES:

- 250 programme furnace
- Motorized pressing system
- Pressing process under continuous control
- 9 programmable pressing levels
- Single or multi-points injection
- Possibility to press in large-sized rings (up to 8 cm of diameter)
- 50 pressing programs
- Available with oil free vacuum pump.

TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz
Dimension	: 26 Kg
Weight	: 400 x 320 x 460 mm



Artis

Furnace for firing high and low fusing ceramics >>

MAIN FEATURES:

- 100 free programs.
- New 'Ugin' Technology System heating chamber.
- Exceptional firing quality.
- Quartz-sheathed resistance.
- Precise control.
- Electrical safeguards.
- Real-time program correction.
- Intuitive navigation (5 languages)
- Excellent readability.
- Elegant and understated design.
- Small Size.
- Wide brushed aluminium shelf.
- Printer connection.
- Available with oil free vacuum pump.

TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz, 1200W
Weight	: 19.5 kg
Dimension	: 520 x 300 x 420 mm

**NEW MODEL
OF LECTRA**



Ceramic Furnace & Vacuum Pump

Ellipse
Ceramic Furnace
MAIN FEATURES:

- 100 free programs
- Heating element protected by a quartz coiled tube UTS.
- Electronic display in 5 languages.
- Controlled cooling.
- Automatic Standby.
- Electronic protection.
- Available with oil free vacuum pump.

TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz, 1300W
Weight	: 19 kg
Dimension	: 480 x 280 x 380 mm


FV Plus
Ceramic Furnace >>

MAIN FEATURES:

- Metal Ceramic Furnace
- Glass Tube Quartz Muffle
- Numeric Display
- 50 Free Programs
- Fast Cooling
- Independent Vacuum & Burning Time
- Energy Saving
- Independent Ascending & Descending Speed
- Compact & Space Saving Design
- Available with oil free vacuum pump.

TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz, 1400W
Weight	: 11 kg
Dimension	: 285 x 320 x 390 mm


AMPRO
Vacuum Pump
Oil-free Piston Pump >>

MAIN FEATURES:

- Outflow : 51 L/min (50 Hz)
55 L/min (60 Hz)
- Limited vacuum : 740 mm/Hg - 29, 13 in/Hg
- Rate current : 1.95 or 2.38 A
- Capacitor : 8 - 30 µf
- Noise Level : 56 dB


TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz, 270W
Weight	: 8.7 kg
Dimension	: 210 x 155 x 315 mm

Burnout Furnaces

Programix 50
Burnout / Pre-heating Furnace
MAIN FEATURES:

- High thermal yield
- Reduced energy consumption
- Compact size
- Modular electronic valve
- Protection against corrosion (electro galvanized steel and epoxy paint)
- Stainless steel door
- Monobloc cast ceramic muffle
- No. of programs 30
- Program regulation 0 to 19 in 3 level stages
- Program regulation 20 to 29 in 6 level stages

TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz, 2000 W
Weight	: 34 kg
Dimension	: 500 x 390 x 440 mm


AMPRO
Burnout Furnace >>

MAIN FEATURES:

Burnout Furnace with Ceramic Fibre insulated and made with thick Mild Steel Sheet. The furnace is insulated with Ceramic Fibre Blanket heavy duty. Heating element made of Kanthal wire spiral type. Heavy duty ceramic fibre blanket is used for better insulation and avoid heat losses. All furnaces are provided with indicator. The furnace is quick heating type and time of heating can be varied as per customers requirement. The apparatus is complete with silver thermal fuse with two meter power lead to operated at 220 / 230 volts.

- Maximum Temperature 1200 °C
- Working Temperature 1150 °C
- With Digital Display Cum Controller

TECHNICAL DATA

Power Supply	: 230 V- 50/60Hz, 2.5 K.W.
Weight	: 30 kg
Dimension	: 100 x 100 x 225 mm






Ugivac
Vacuum Mixer

MAIN FEATURES:

- Electronic Timer
- Digital Display
- Adjustable Vacuum Level
- Automatic coupling system
- Rotation speed 400 tr/min
- Delivered with 2 bowls : 250 cc & 500 cc
- In-built filter

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 100 W
Weight	: 11 kg
Dimension	: 180 x 255 x 230 mm



Mix-R
Vacuum Mixer

MAIN FEATURES:

- Wall mounted vacuum mixer.
- Adjustable mixing speed and possibility to reverse the paddle direction providing smoother and more homogeneous mixing.
- Programmable timer with buzzer at the end of the cycle.
- Accurate and constant vacuum mixing
- Supplied with 500 cc bowl.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 140W
Weight	: 8.7 kg
Dimension	: 220 x 250 x 430 mm



MAIN FEATURES:

- Bubble-free mixing of investments and stones under vacuum.
- Two mixing pot sizes available for investments, silicones and dental stones.
- Digital timer for setting mixing time provided.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz
--------------	------------------


W.A.F.I.S
Water Air Filtering System

Dust particles will deposit inside the machine, not in your lungs!

The WAFIS system is proposed in different versions to suit with any sandblaster on the market thus updating your appliance according to the Regulations in terms of environmental hygiene.

THE ADVANTAGES:

- Complete elimination of dust:
99,8% assured
- No maintenance cost:
water only to be replaced
- Ridiculous running cost:
very low air consumption
- Excellent quality/price ratio:
no alternative system is less in costs
- Full compliance with regulations in force:
respect for environment hygiene


Base 3 Sand Blaster
With two tanks with recycling nozzle
MAIN FEATURES:

- Divesting Nozzle for recyclable abrasive 3.0 mm nozzle
- EASY micro blasting tanks for disposable media fitted with 0.8 mm nozzle (0.5 - 0.8 - 1.2 - 2.0 mm nozzles available on request)
- Air Filter
- Pressure gauge and foot control
- Pre-fitted for connection to WAFIS and PRO-3. Suction system

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 60 W
Weight	: 11.8 kg
Dimension	: 415 x 400 x 445 mm


AMPRO
Sand Blaster

MAIN FEATURES:

- Sand Blaster (Dual tank) for abrasion
- Free working on casting with nozzle allowing contamination free smooth finish of articles,
- Built with light-weight material,
- Exhaust filter for pollution free atmosphere and large glass window for easy viewing.
- Supplied with 2 nozzles

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz
Weight	: 9.5 kg
Dimension	: 375 x 400 x 425 mm

Suctions for Sandblaster

Pro 3 Shake
Suction for Sandblaster
MAIN FEATURES:

- High performance electrical vacuum cleaner
- Supporting 1 or 2 units.
- Power Supply is ensured also directly on the bench working areas.
- Electromagnetic stirrer SHAKE.
- An unique application which will repeatedly clean the filter reducing the need to access the inside of the unit.
- The 3-stage-filtration concept ensures a perfect protection of the engine.
- Paper bags are sold in packs of 10.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 1200W
Weight	: 13.3 kg
Dimension	: 180 x 255 x 230 mm


Pro 3
Suction for Sand Blaster
MAIN FEATURES:

- More Compact
- Suitable for reasonable suction volumes. It can also be connected to two units at the same time. Paper bags are sold in packs of 10.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 1000W
Weight	: 9.8 kg
Dimension	: 380 x 380 x 620

Surveyor / Milling

ISO Color Top
Surveyor & Milling Unit
MAIN FEATURES:

- With air
- With light
- With lifting device
- Double-joint arm to support any micro motor
- Turbine or hot spatula
- Most suitable to carry out milling operations
- Complete with model-holder base
- Stainless steel working plan.
- Also available without air, light & lifting device

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 6 w Revolving halogen lamp, 6W absorption
Weight	: 7 kg
Dimension	: 200 x 340 x 330 mm


Surveyors / Parallelometers
Surveyors "B" Type
MAIN FEATURES:

- 1 double-joint arm with pantograph movement ideal for the model planning stage and layout of reference points

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz
Weight	: 4 kg
Dimension	: 190 x 210 x 310 mm





<< **Model Trimmer Dry**
(Single Disk)


MAIN FEATURES:

- Waterless
- Speed 1400 rpm.
- Electric socket for automatic connection to the extraction system.
- Long life diamond coated disc.
- Average-high speed assuring a good polishing effect of the surface whilst cutting down the possibility for the abrasive disc to gum.
- Electromagnetic safety switch avoiding the machine from starting up accidentally after a voltage drop.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 600W
Weight	: 15.7 kg
Dimension	: 320 x 430 x 350 mm



>> **Model Trimmer Wet**
(Single Disk)

MAIN FEATURES:

- Wet model trimmer
- Speed 1400rpm.
- Load solenoid valve adjusting the quantity of water.
- Bakelite disc is reversible.
- A diamond coated disc can be installed at any time.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 600W
Weight	: 15.7 kg
Dimension	: 320 x 430 x 350 mm



>> **AMPRO**
Model Trimmer Wet
(Single Disk)

Model Trimmer with 10" Coarse Cutting imported stone wheel.

TECHNICAL DATA

Power Supply	: 230V AC, 50/60 Hz, 600W
Weight	: 15.7 kg
Dimension	: 320 x 430 x 350 mm



»»

AMPRO

Model Trimmer Wet (Double Disk)

MAIN FEATURES:

- 2800 RPM Double Motors 1 HP each
- Heavy duty equipped with two ball bearing each motor
- Single speed,
- Continuous rating with fixed on e imported 10" coarse cutting stone and one 10" Indian polishing stone wheel
- Heavy base
- Inlet-outlet water flow system with PVC pipe 2 mtrs. Each

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz



»»

GAMBERINI srl

Italy

Model Trimmer Dry (Double Disk)

MAIN FEATURES:

- 2 independent ventilated motors (450 Watts/1400 Rpm).
- Complete with solenoid valve
- Fine water regulation and orthodontic accessories

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz, 2 x 450W
 Weight : 37 kg
 Dimension : 530 x 420 x 400 mm



««

AMPRO

Palatial / Arch Trimmer

MAIN FEATURES:

- Trims the inside and outside surfaces of plaster models with carbide or diamond tools mounted in the spindle that rotates at 5,000 RPM.
- Cuts rapidly and leaves a smooth finish.

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz



>>

AMPRO High Speed Alloy Grinder

MAIN FEATURES:

- Eagle Brand R2 high speed drilling bit and imported motor.
- All kinds of cutting and polishing for various pint sized metal, resin, or plastic work pieces and craftwork.
- Integrated drilling bit and motor at the same time
- Quickly replace various cutting pins, bits and cutting pieces.
- High rotation speed, fast cutting and grinding and low noise.
- Ideal specialized device for dental technicians.

TECHNICAL DATA

Power Supply	: 230V AC, 50/60 Hz, 250W
Supply Voltage	: 220V, 50Hz / 110V 60Hz
Rotation Speed	: 17500 r.p.m
Weight	: 18 kg
Dimension	: 415x380x330 mm



>>

AMPRO Dental Lathe

MAIN FEATURES:

- ½ HP 2800 R.P.M heavy duty rating
- 3 attachment (Buffer, Stone and Bur) & D.B.B Cable Arm.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz
--------------	------------------



>>

GAMBERINI srl Italy Aspiration Hood Polishing Lathe with suction

MAIN FEATURES:

- Filter tank
- Halogen lighting
- Stainless steel work surface
- Rubber splash bins

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, Aspirator: 180W, vLathe: 300W-1400/2800 Rpm
Weight	: 46 kg
Dimension	: 670 x 500 x 515 mm



>>

AMPRO Dental Hanging Engine

MAIN FEATURES:

- 18,000 R.P.M
- Imported three B.B cable arm
- 5 speed foot controller.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 400W
--------------	------------------------

Laser Pin / Die Cutting


Laser Pin 
Pindex Machine Laser Pin-Drilling Unit
MAIN FEATURES:

- Precision laser pinhole drilling unit.
- The accurate drill hole positioning system and the easily adjusting working position offer the best warranty for excellent results.
- Zeiser type drill supplied with.

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz, 150W
 Weight : 15 kg
 Dimension : 250 x 260 x 420 mm


AMPRO 
Laser Pin
Pindex Machine Laser Pin-Drilling Unit
MAIN FEATURES:

- Laser Pin for super-precision laser adjustment point
- Adjustable diameter (between 0.5 and 3 mm)
- Clear view of the model.

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz, 150W
 Weight : 15 Kg
 Dimension : 250 x 260 x 420 mm

GAMBERINI srl
 Italy

Die Cutting Unit with Suction 
Trunk Separator
MAIN FEATURES:

- Diamond blade (Ø80 x Ø16 x 0.2mm)
- Internal suction system
- Halogen lighting
- Dust filter
- Standard model support
- Internal Suction
- 4800 rpm saw motor
- Inbuilt electromagnet
- Dual security

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz, 80W
 Weight : 22.2 kg
 Dimension : 260 x 400 x 365mm





Pressojet Steam Cleaner

MAIN FEATURES:

- Water level safety device.
- Audible signal.
- indicator light and electrical cut out.
- Steam gun outlet and additional dry vapor outlet via fixed nozzle
- Efficiently cleans tour denture work in complete safety.
- Stainless steel water tank of 4 liters (Operating capacity 3 liters)
- Pressure 3.5 bars.
- Foot control for electro-valve.
- Water level indicator.

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz, 1600W
 Weight : 12 Kgs
 Dimension : 320 x 320 x 330 mm



Ducagel - ZI Gel Duplicating Machine

- Heating Homogeneity, Eliminating Overheating Risks
- 2 Possible Compounding Procedures-Continuous/Intermittent
- 10 Programs

TECHNICAL DATA

Power Supply : 230V 50/60 Hz, 800W
 Rotation Speed : 30 tr/min
 Capacity : 6 kg
 Weight : 20 kg
 Dimension : 325 x 460 x 410 mm



Hydraulic Flask Press


MAIN FEATURES:

- Described as reliable.
- Easy to use with sturdy
- Durable construction.
- Can accommodate up to three flasks
- Pressure of 400bar.
- The oleo dynamic circuit is absolutely leak-proof and no maintenance is required.

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz
 Weight : 22 kg
 Dimension : 200 x 200 x 550 mm



AMPRO Electrolytic Polishing Unit

MAIN FEATURES:

- Applications for non oxidizing dental alloys.
- Provides surfaces with bright, glossy and durable luster.

TECHNICAL DATA

Power Supply	: 220 V 50/60 Hz
Weight	: 9 Kgs
Dimension	: 175 x 245 x 150 mm



AMPRO Magnifying Lamp Bench Lighting

MAIN FEATURES:

- Lamp fitted with large
- Power Supply full (up to 3 magnifying options) lens allowing to be extended up to 180mm.
- Two 9W lamps ensuring 6400°K, thus offering white natural light.
- Supplied complete with hangers to fix the lamp either on the bench or to the wall.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 2 x 9W
Weight	: 5.5 kg
Dimension	: 300 x 200 x 50 mm



AMPRO Vibrator

MAIN FEATURES:

- Vibrator capacity 3-4 Casting Cylinder,
- Removable rubber plate for easy cleaning,
- Fitted with controller for variable vibrational frequencies.

TECHNICAL DATA

Power Supply	: 230V AC, 50/60 Hz
--------------	---------------------



AMPRO Plaster Dispenser

MAIN FEATURES:

- Electric controlled unit for dispensing required quality of plaster at steady rate with On/Off control switch.
- Capacity-10Kg. Can be wall mounted

TECHNICAL DATA

Power Supply	: 230 V AC, 50/60 Hz
--------------	----------------------

AMPRO Phantom Head Simulator Unit >>

MAIN FEATURES:

- Phantom head Simulator Unit includes
- Motorized Phantom Head with in-built Halogen light.
- Airtor control box with controls for Airtor.
- Micromotor (with 2 hand pieces) and 3-way syringe.
- Pneumatic Stool, Phantom head with half fibre glass body.
- Motorized up & down movement of body.
- Touch Pad control panel.
- Imported glass reflector.
- Typodont Jaw Set with 32 ivory teeth.
- Enamel Finish.

TECHNICAL DATA

Power Supply : 230 V 50/60 Hz



>> AMPRO Phantom Head Table Model

MAIN FEATURES:

- 32 Typodont Ivory Teeth
- Face Mask With Drainage System
- Metal Articulator
- Adjustable Arm Height For Table
- All Natural Movement Posterior / Lateral Anterior



>> AMPRO Jaw Set With Plastic Articulator

MAIN FEATURES:

- 32 Typodont Ivory Teeth
- Plastic Articulator
- Soft Gums.



>> AMPRO Jaw Set With Metal Articulator

MAIN FEATURES:

- 32 Typodont Ivory Teeth
- Metal Articulator
- Soft Gums.



Linea 1 Workstation For Clinic

Composition	W2 + W4 + L(55) + W20/W21
Accessories	LAV.CER + GM/S
W2	Seven Drawers Model
W4	Sink model W4
L(55)	
W21	

Modular unit on wheels specially designed to contain the instruments, consumables, trays, towels etc. that are required for the dental chair.

Without top

W20	Modular open unit for autoclave. Without top
LAV.CER	Ceramic basin over-top
GM/S	Mixer set for cold-warm water
PM/S	White laminated plastic top -30 mm thick

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz
Dimension	: 208 x 53.5 x H86.5 cm



Sole 1 Workstation For Lab

MAIN FEATURES:

- ASOLE1 Workstation for one person.
- Worktop made of laminated plastic.
- Suction Set Complete with variable speed control
- Stainless Steel worktop protection
- Adjustable electronic left hand lamp, 55w Power Supply, 98 Ra, 5300 K. Two drawers complete with internal plastic trays. Central drawer unit equipped with three drawers of which one is lockable. Gas tap and pipe work, Airgun & Micro Motor support.
- Four electrical sockets
- The entire structure, working positions and drawer unit are manufactured in galvanized steel in accordance with standard EN-10152, a first choice for its quality and uniform grade. Painting is completed with enamel powder, polymerized in oven at 180 C, and provides a flexible, strong and homogeneous surface resistant to abrasion, acids, alkalis and solvents.

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz, 55W
Dimension	: 110.5 x 60 cm



Linea 3 Workstation For Clinic

Composition	W19 + W13 + W1 + W20
Accessories	LAV.CER+GM/S
W19	Modular unit with door and drawer. The unit includes a shelf without top
W13	Corner unit with door, and has provision for a wash basin. Without top
W1	
W20	Modular open unit for autoclave. Without top
LAV.CER	Ceramic basin over-top
GM/S	Mixer set for cold-warm water
PM/S	White laminated plastic top-30 mm thick

TECHNICAL DATA

Power Supply	: 230 V 50/60 Hz
Dimension	: 121 x 171 x 53.5 x H86.5 cm


Anand®
EUROVAC Suction Unit
Portable (Single jar)
MAIN FEATURES:

- Anand's Real Portable Unit with Shock Resistant ABS Case.
- Vacuum: -550mm Hg \pm 10 at 22 ltrs / min,
- Oil-free diaphragm pump. Autoclavable 2 ltr PC jar with over flow safety jar.

TECHNICAL DATA

Power Supply	: 220V AC, 50/60 Hz, 125 watts (110V on request)
Weight	: 7.5 kg
Dimension	: 20 x 42 x 36 cms

Anand®
MB-36 Suction Unit >>
Heavy Duty (Double jar)

- Electric Operated Suction Unit. Spot Welded, Oven baked cabinet.
- Vacuum: -700 mm Hg \pm 10 regulable, 25 Ltrs, rotary vane pump
- Jars : Wide mouthed 2 x 2 Ltrs PC jars with over flow safety device.

TECHNICAL DATA

Power Supply	: 230V, 50 Hz, 2 \pm 0.5 Amps, 200 watts (110V on request)
Weight	: 27 kg
Dimension	: 50 x 38 x 70cms





WOL-CERAM[®]
 WOL-CERAM ONE[®]
 Electroplating CNC Centre

MAIN FEATURES:

- Guarantees your independence, remaining competitive working economically and with precision
- Equipped with 2 axis, touch screen control element
- High performance ceramic frameworks made of alumina and zirconia
- 9 resident programs for bridges and 1 resident program for implant framework
- Gold, Silver, non-ferrous metal alloys
- SPS-controlling
- Productivity 5 units per hour.


TECHNICAL DATA

Power Supply : 230V AC, 50/60 Hz

Wol-DENT Gmbh guarantees the possibility to carry out both metal and ceramic bridges by electrophoresis in a quick and precise way.



Wol-ceram ELC/ELM allows the production of both metal-free and PFM alloy structures.



Wol-ceram guarantees the high quality of metal powder (obtained by atomising), produced thanks to its decade-long experience with the most modern production system. Certified quality



Wol-ceram means: reduced manual work, no casting bottoms or waste of material, repair and modification, always new material (important when using alloys: no recasting of metal), unbelievable fit, total autonomy of the lab (not depending by outsourced work).



The powders are the basic component of the material to be used with the electrophoresis. The quality supplied by Wol-ceram together with the patented working technique by WOL-DENT assures a high technology execution level of an electrophoresis cycle which leads to the production of metal crowns and alumina bridges in just 5 minute time.

WOL-CERAM[®]
 WOL-TEMP 1200 M
 Sintering Furnace



Consumables


Performance Plus Press is a new type of low fusing; leucite reinforced pressable dental porcelain that allows both staining and layering technique. Performance Plus Press can be pressed directly on a metal structure made in SINTERCAST (single crowns), in precious and non precious alloys (single crowns and bridges).



Ugicast 24 for precision investment of precious and non-precious alloys. Main features include, plaster less, controlled expansion, highly precise casting, smooth surfaces and fine granulometry. (Box of 6.08 kgs) 38 pouches of 60gms each with 1ltr liquid.



Ugirex III Non precious alloy for ceramic (without beryllium). Pack of 100gms. Also available - Alloys for Cr-Co and Ni-Cr alloy for crown and bridge.



Uginvest 38 is an investment for Chrome - Cobalt with smooth vitrified surface & High precision of casting. Available in white color.

Expansion liquid Uginvest plus. (Box of 15kgs) 75 pouches of 200gms each with 1ltr liquid.



Crucible CARBON for Ducatron Serie 3. (Sold by 3 units)/Kerr type (Sold by 5 units)



Ring liner for expansion (without asbestos) - In the form of strips. - Length 25 m, height 50 mm, thickness 1 mm




Alumina Sand 50/100/150/250 microns. Its Utilization Surface retention is on all types of Air pressure ¾ bars Diam, Nozzle 0.8mm to 2mm. Available in 10 Kg packs



Wax/Ugifil for precision investment of precious and non-precious alloys. Main features include, plaster less, controlled expansion, highly precise casting, smooth surfaces and fine granulometry (Box of 453gms). Also available in wire form (Box of 250gms).

AMPRO undertakes Turnkey Set up of Dental Colleges,
Dental Laboratories and Dental Clinics.
Our Engineers and Consultants are well equipped to design
as per your requirements and space availability.



SOME OF OUR ESTEEMED CLIENTS

ALL INDIA INSTT. OF MEDICAL SCIENCE	NEW DELHI
SIR GANGA RAM HOSPITAL	NEW DELHI
ITS CENTRE FOR DENTAL STUDIES & RESEARCH	GAZIABAD, U.P.
P.D. MEMORIAL RELIGIOUS & EDUCATIONAL ASSOCIATION (REGD.)	JHAJJAR, HARYANA
MANAV RACHNA EDUCATION SOCIETY.	FARIDABAD. HARYANA
D.J. COLLEGE OF DENTAL SCIENCES	GAZIABAD, U.P.
TEERANKER MAHAVIR DENTAL COLLEGE & HOSPITAL	MORADABAD, U.P.
SEEMA DENTAL COLLEGE & HOSPITAL	RISHIKESH
KOTHIWAL DENTAL COLLEGE & RESEARCH CENTRE	MORADABAD, U.P.
S.R.M. DENTAL COLLEGE	CHENNAI, TAMIL NADU
MAHATMA GANDHI NATIONAL INSTT. OF MEDICAL COLLEGE	TONK ROAD, JAIPUR
AL - AMEEN DENTAL COLLEGE	BIJAPUR, KARNATAKA
MILITRY DENTAL COLLEGE, ARMED FORCES	PATHANKOT
SARASWATI DENTAL COLLEGE	LUCKNOW
K.M. SHAH DENTAL COLLEGE & HOSPITAL	VADODARA
JAMIA MILIA ISLAMIYA DENTAL COLLEGE	NEW DELHI
SARDAR PATEL DENTAL COLLEGE	LUCKNOW
MODERN DENTAL COLLEGE	INDORE
K.D. DENTAL COLLEGE	MATHURA
K.G.M.C. COLLEGE	LUCKNOW
PEOPLE'S COLLEGE OF DENTAL SCIENCES & RESEARCH CENTRE	BHOPAL
GYAN SAGAR DENTAL COLLEGE	PUNJAB
P.M. NADAGOUDA MEMORIAL DENTAL COLLEGE & HOSPITAL	BAGALKOT, KARNATKA
GINNI MODI DENTAL COLLEGE	MODINAGAR
SANTOSH DENTAL COLLEGE	GHAZIABAD, U.P.
HARSARAN DASS DENTAL COLLEGE	GHAZIABAD, U.P.
SRI GOVIND TRICENTENARY DENTAL COLLEGE HOSPITAL & RESEARCH INSTITURE	GURGAON
GIAN SAGAR DENTAL COLLEGE	CHANDIGARH
RUNGTA DENTAL COLLEGE	KOHOKA, BHILAI
SWAMI DEVI DYAL DENTAL COLLEGE	BARWALA, PANCHKULA
COMMAND MILITARY DENTAL COLLEGE	PUNE, MAHARASHTRA
COMMAND MILITARY DENTAL COLLEGE	CHANDIGARH
MILITARY DENTAL COLLEGE	ANDHRA PRADESH
GOVT. DENTAL COLLEGE	GUWAHATI
GOVT. DENTAL COLLEGE	RAIPUR
HIMACHAL DENTAL COLLEGE	HIMACHAL PRADESH
HITKARNI DENTAL COLLEGE & HOSPITAL	JABALPUR, M.P.
BANKEY BIHARI DENTAL COLLGEG & RESEARCH CENTRE	GAZIABAD, U.P
ARODENT DENTAL LAB	NEW DELHI
WESTMAN INTERNATIONAL DENTAL LAB	NEW DELHI
GANI DENTAL LAB	CHENNAI
SRI SAI DENTAL LAB	HYDERABAD
M.DENT DENTAL LABORATOR	AHEMDABAD
KATARA DENTAL LAB	PUNE
DR. P. RAJKUMAR	TAMIL NADU

And many more...