

SAIKANG[®]

LEADING BRAND OF CHINA
MEDICAL FURNITURE

JIANGSU SAIKANG MEDICAL EQUIPMENT CO., LTD.

www.saikangmedical.com

Stock Code: 870098



SAIKANG®

LEADING BRAND OF CHINA
MEDICAL FURNITURE



LEADING BRAND OF CHINA
MEDICAL FURNITURE

RACING AGAINST TIME
HEALTHIER YOUR LIFE



CONTENTS

01	Strength of Saikang	09	Hospital Bed	69	Medical Bridge, Ceiling Pendant
03	Technology of Saikang	31	Homecare Bed	74	Mattress
05	Innovation of Saikang	33	Children Bed	76	Bedside Table
07	Quality of Saikang	39	Baby Crib	77	Overbed Table
		41	Spinal Therapy Bed	78	Medical Chair
		43	Examination Table	84	I.V. Stand
		47	Obstetric Bed	85	Stretcher Trolley
		52	Gynecological Table	86	Patient Trolley
		53	Operating Table	91	Medical Trolley
		65	LED Operation Lamp		

STRENGTH OF SAIKANG

24

24 years focusing on design, production, service and technical consultation of medical care equipment

150

Located in Zhangjiagang City, 150km nearby Shanghai Economic Center

10000000

Saikang delivered 1,000,000 sets of products globally. Products all over Asia, Middle East, North America and Africa

400

There are more than 400 employees worldwide, more than 50 qualified engineers, and more than 100 trained technicians

245

Equiped with 245 sets of various automatic production facilities, Saikang follow strict standard control system on every production process, enable the traceable of each part

50000

The finished phase I production area covers more than 30,000 square meters, and phase II will be more than 20,000 square meters of advanced production area

31

Saikang offer comprehensive medical solutions with 31 product categories

562

Safe, efficient and easy-to-use medical solution including emergency, nursing and ICU wards of more than 500 series of products

32

With CE certification, ISO9001/13485 quality system certification, U.S. FDA and other 32 certifications

189

Ownership of 189 patents, continuously leading new standards in defining product quality and service efficiency

150

Since the establishment of the company, products and solutions have been used in more than 150 countries

1000

Saikang products are committed to creating a modern nursing environment, accumulatively serving more than 1,000 medical and health care organizations



LEADING BRAND OF MEDICAL FURNITURE

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in US, Dubai, Russia, Brazil, Thailand and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 149 countries and regions in the world.

TECHNOLOGY OF SAIKANG

SAIKANG insists on providing reliable quality, reshapes the "Made in China" products in worldwide market. SAIKANG has obtained ISO9001/13485 quality system certification, CE, and approved by the U.S. Food & Drug Administration (FDA) for its products. SAIKANG has also obtained EU safety certifications, and US FDA 510(k) clearance and Russian FDA market access license.

Alibaba.com Assessed Supplier

Supplier Assessment
Certificate Number: 779799 P-F
Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司
For details, please refer to <http://saikang.alibaba.com>

CE Technical Documentation Review Report

Manufacturer: Jiangsu Saikang Medical Equipment Co., Ltd.
Report Number: 1508465791
Examination item: Operating Tables
Product(s): SPT series, A series, SPT series
Type/Model: Class 1, Class 2
Classification: Class 1, Class 2
Review Result: Passed

This is to certify that the Quality Management System of Jiangsu Saikang Medical Equipment Co., Ltd. has been assessed and approved by the Quality Management System of the International Organization for Standardization (ISO) and the Quality Management System of the International Organization for Standardization (ISO).





INNOVATION OF SAIKANG

Design Innovation



SAIKANG medical designers have always maintain close communication with medical service providers, to help them providing better solutions to their patients. With long-standing experience in product structure and functional experience, SAIKANG continue to forge ahead with innovation and play a vital role in helping medical service providers to improve the quality of care.

Research Innovation



SAIKANG medical owns independent research and development center. Since its inception, every year SAIKANG invested 10% of its annual sales in R&D, launched 10 new products, and applied at least 10 patented technologies on an average basis, which lead SAIKANG to continuously define the quality of medical products and new standards of efficiency.

Manufacturing Innovation



SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality.

Maintenance Innovation



SAIKANG medical have more than 50 engineers, and qualified service term that composed of professional subcontractors, trained technicians, and after-sale staff. They are experienced and have excellent skills , not only to respond each service request on time, but also precaution every detail in equipment operation , always focus on products safely.

Service Innovation



SAIKANG medical have thorough service architecture: 7days x 24hours online customer call center, complete spare parts warehouse, advanced CRM system, precautionary service, etc. No matter where customer is, SAIKANG service will be always close at hand.

 Doctor-patient communication system

Rehabilitation guidance, follow-up, tracking, reminders, allowing patients to communicate with doctors through text, voice, and video anytime, anywhere, creating a mobile medical platform for in-depth doctor-patient communication.



 Real-time vital sign monitoring

Monitor heart rate, blood pressure, blood oxygen, respiration, temperature, pulse, etc.



Remote Monitoring System

Real-time video monitoring and audio-video conversation.




**Vital Sign Monitoring
Patient 360 Monitoring
Infusion Monitoring
Video Monitoring**




Equipped with Smart Hospital Bed System

Clinical data collection, screening, integration. Generate medical big data to assist decision-making.



Control system





75°±10°
Back-rest

35°±10°
Knee-rest

400mm
Height
Adjustment
530-930mm

15°±2°
Reverse
Trendelenburg

15°±2°
Trendelenburg

17°±2°
Lateral Tilt
Left

17°±2°
Lateral Tilt
Right

weight
Weighing System

CPR
Instant Release

Cardiac
Chair Position

0-90°
Angle Indicator

75+15mm
Auto-regression

30°
Position Line

LINAK
Motor

Brake Indicator

Rail Lock
Indicator

Z7z
ICU Smart Bed
SKZH-ABCZ

Technical Parameters

External Length	2127mm
External Width	1070mm
Bed Height Min.	530mm
Bed Height Max.	930mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ LINAK Motor	2pcs
✓ Telescopic Column Motor	3pcs
✓ 6" Double Side Secure® Castors, Central Locking	4pcs
✓ PP Side Rail with Integrated LCD Control Panel	1set
✓ PP Bed Ends with Safe Lock	1set
✓ X-ray Transparent Phenolic Resin Bed Platform	1set
✓ Patent Angle Indicator	4pcs
✓ Φ85mm Wall Bumper Wheel	4pcs
✓ Trapeze Socket	2pcs
✓ IV Pole	1set
✓ IV Pole Prevision	8pcs
✓ Drainage Hook	4pcs
✓ Central Locking Pedal	2pcs
✓ PP Mattress Retainer	6pcs
✓ Nurse Cotroller	1pc
✓ Emergency Stop Button	1pc
✓ Night Light	1set
✓ Manual CPR	
✓ Electric CPR	
✓ Weighing System	



70°±10°
Back-rest

28°±10°
Knee-rest

450mm
Height
Adjustment
400-850mm

17°±2°
Reverse
Trendelenburg

17°±2°
Trendelenburg

Cardiac
Chair Position

200mm
Bed Extension

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

67±57mm
Auto-regression

30°
Position Line

Built-in LED
Torch

LINAK
Motor

Nurse Call

Siderail
Unlocked

Brake
Unlocked

Lowest
Position

Out-of-bed
Position Warnings

C8c
Electric ICU Bed
SKD-CV

Technical Parameters

External Length	2180mm
External Width	995mm
Bed Height Min.	400mm
Bed Height Max.	850mm
Mattress Platform	1995×860mm
Safe Working Load	250kg

Technical Configuration

✓ LINAK Motor	4pcs	✓ Nurse Cotroller	1pc
✓ 5" Double Side Castors with directional wheel function	4pcs	✓ Cardiac Chair Position	
✓ PP Removable C-series Headboard and Footboard Φ113 Bumper Wheels	1set	✓ Electric CPR	
✓ PP 4-pcs C-Series Siderails	1set	✓ Manual CPR	
✓ Head-end Siderails with Integrated Controller and Warning Lights	1set	✓ Auto Break at 30°	
✓ Foot-end Siderails with Glasses Holder, LCD Screen and 15W Wireless Charger	1set	✓ Emergency Stop Button	
✓ PP Detachable Bed Platform With Integrated Mattress Retainer	1set	✓ Night Light	
✓ Mechanical Lever on Leg Board for Calf Raise	1set	✓ Bed Extension	
✓ Backrest with X-ray Cassette Holder	1pc	✓ Weighing System	
✓ Patent Angle Indicator for Back-rest and Trendelenburg Position	4pcs	✓ Bed End with Projection Lights	
✓ IV Pole 4 Hooks	1set	✓ Nurse Call	
✓ IV Pole Prevision	4pcs	✓ Siderail Unlocked	
✓ Drainage Hook	2pcs	✓ Brake Unlocked	
✓ Φ32mm Monkey Pole Holes	2pcs	✓ Lowest Position	
✓ Hand Cotroller	1pc	✓ Out-of-bed Position Warnings	



75°±10°
Back-rest

20°±10°
Knee-rest

400mm
Height
Adjustment
550-950mm

18°±2°
Reverse
Trendelenburg

18°±2°
Trendelenburg

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

75mm
Auto-regression

200mm
Bed Extension

Cardiac
Chair Position

LINAK
Motor

X9x
Electric ICU Bed
SKD-CV

Technical Parameters

External Length	2120-2320mm
External Width	1010mm
Bed Height Min.	550mm
Bed Height Max.	950mm
Mattress Platform	1890-2090×820mm
Safe Working Load	250kg

Optional

☐ Monkey Pole

Technical Configuration

✓ LINAK Motor	6pcs
✓ 6" Double Side Castors, Central Locking	4pcs
✓ PP Side Rail with Integrated Control Panel	1set
✓ PP Bed Ends One Step Removable	1set
✓ Ø85mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	4pcs
✓ 5th Directional Wheel Pedal	2pcs
✓ Central Locking Pedal	2pcs
✓ X-ray Translucent Mattress Platform	1set
✓ Back-up Battery	1set
✓ Night Light	1set
✓ Bed Push Handle for Emergency Transfer	1set
✓ Manual CPR	
✓ Electric CPR	
✓ Weighing System With LCD Screen	



75°±10°
Back-rest

33°±10°
Knee-rest

450mm
Height
Adjustment
390-840mm

17°±2°
Reverse
Trendelenburg

17°±2°
Trendelenburg

Cardiac
Chair Position

200mm
Bed Extension

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

70+40mm
Auto-regression

30°
Position Line

Built-in LED
Torch

LINAK
Motor

Oasis Series

Y8y
Electric ICU Bed
SKD-CV

Technical Parameters

External Length	2150mm
External Width	1060mm
Bed Height Min.	390mm
Bed Height Max.	840mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Optional

- ☐ Monkey Pole
- ☐ Oxygen Cylinder Holder
- ☐ Basket

Advantages

- ★ Double Auto-regression
- ★ Super-low Bed Position 390mm
- ★ PP Bed Platform, Removable & Washable
- ★ Inserted-type PP Head Foot Board, One Step Release With Safe Lock, Patent Design
- ★ PP Safety Guardrail With Two Steps Unlocking System, Patent Design
- ★ Back Rest 30degree Position, Reduce The Chance Of Lung Infection
- ★ Built-in LED Torch

Technical Configuration

✓ LINAK Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ PP Oasis Side Rail Extra Height	1set
✓ PP Oasis Bed Ends with Safe Lock	1set
✓ PP Bed Platform, Removable & Washable	5pcs
✓ X-ray Transparent Phenolic Resin Back Platform	1set
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø113mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	6pcs
✓ Three-grade Central Locking Pedal	2pcs
✓ ABS Mattress Retainer	2pcs
✓ Manual CPR	
✓ Electric CPR	
✓ Nurse Controller	
✓ Hand Controller	
✓ Weighing System	
✓ Night Light	
✓ Bed Extension 200mm	
✓ Trapeze Socket	
✓ Back Up Battery	
✓ Fif th and Six th Siderails	



75°±10°
Back-rest

35°±10°
Knee-rest

290mm
Height
Adjustment
480-770mm

12°±2°
Reverse
Trendelenburg

12°±2°
Trendelenburg

Cardiac
Chair Position

CPR
Instant Release

0-90°
Angle Indicator

30°
Position Line

LINAK
Motor

Oasis Series

Y8t
Electric ICU Bed
SKD-C

Technical Parameters

External Length	2140mm
External Width	1060mm
Bed Height Min.	480mm
Bed Height Max.	770mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ LINAK Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ PP Oasis Side Rail Extra Height	1set
✓ PP Oasis Bed Ends with Safe Lock	1set
✓ Stamping Forming Bed Platform	1pc
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø113mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ Emergency Stop Button	1pc
✓ Manual CPR	
✓ Electric CPR	
✓ Lowest position indicator	

Optional

☐ Hand Controller
☐ Weighing System
☐ Nurse Control Panel
☐ Cardiac Chair Position
☐ Bed Extension 20cm
☐ Back-up Battery
☐ Directional Wheels Function
☐ ABS Mattress Retainer
☐ Monkey Pole
☐ Balkan Frame
☐ Linen Holder



360° View



75°±10°
Back-rest

35°±10°
Knee-rest

290mm
Height
Adjustment
435-725mm

12°±2°
Reverse
Trendelenburg

12°±2°
Trendelenburg

CPR
Instant Release

0-90°
Angle Indicator

200mm
Bed Extension

LINAK
Motor

Grace Care

D8d-CPR

Electric Bed
SKD-C



Technical Parameters

External Length	2280mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Optional

- ☐ Monkey Pole
- ☐ Oxygen Cylinder Holder
- ☐ Basket

Technical Configuration

☒ LINAK Motor	4pcs
☒ 5" Double Side Castors, Central Locking	4pcs
☒ PP Grace Side Rail with Integrated Control Panel	1 set
☒ PP Grace Bed Ends with Safe Lock	1 set
☒ Intergrated Control Panel (3.5 inch LCD)	1 pc
☒ Phenolic Resin Backrest platform	1 pc
☒ Trapeze Socket	2pcs
☒ IV Pole 4 Hooks	1 set
☒ IV Pole Prevision	4pcs
☒ Wall Bumper Wheel Φ85mm	4pcs
☒ Three-grade Central Locking Pedal	2pcs
☒ Drainage Hook	2pcs
☒ Bed extension 200mm	
☒ Manual CPR	
☒ Electric CPR	
☒ Emergency Stop Button	
☒ Cassette Holder 448×400mm	
☒ Night Light	
☒ Back-up Battery	



75°±10°
Back-rest

28°±10°
Knee-rest

450mm
Height
Adjustment
420-870mm

17°±2°
Reverse
Trendelenburg

17°±2°
Trendelenburg

Cardiac
Chair Position

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

67±47mm
Auto-regression

30°
Position Line

LINAK
Motor

V8v
Electric Bed
SKD-CV

- RAL 210 710 30
- RAL 5023
- RAL 6019
- RAL 3015
- Pale Oak
- White Oak
- Red Oak

Technical Parameters

External Length	2085mm
External Width	1025mm
Bed Height Min.	420mm
Bed Height Max.	870mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Optional

- ☐ Monkey Pole
- ☐ Oxygen Cylinder Holder
- ☐ Basket
- ☐ Mechanical Lever

Technical Configuration

✓ LINAK Motor	4pcs	✓ Nurse Cotroller	1pc
✓ 6" Double Side Castors, Central Locking	4pcs	✓ Trapeze Socket	2pcs
✓ PP Side Rail Extra Height	1 set	✓ Cardiac Chair Position	
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1 set	✓ Electric CPR	
✓ PP Bed Platform with Integrated Mattress Retainer	1 set	✓ Manual CPR	
✓ X-ray Translucent Mattress Platform (Back Section)	1pc	✓ Auto Break at 30°	
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs	✓ Emergency Stop Button	
✓ IV Pole 4 Hooks	1 set	✓ Night Light	
✓ IV Pole Prevision	4pcs	✓ Cassette Holder (Backrest)	
✓ Drainage Hook	2pcs	✓ Bed Extension	
✓ Central Locking Pedal	2pcs	✓ Weighing system	
✓ Hand Cotroller	1pc		

V6v
Electric Bed
SKD-D

75°±10°
Back-rest

46°±10°
Knee-rest

310mm
Height
Adjustment
430-740mm

0-90°
Angle Indicator



Technical Parameters

External Length	2085mm
External Width	1025mm
Bed Height Min.	430mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	3pcs
✓ Hand Controller	1pc
✓ 5" Hardshell Castors	4pcs
✓ PP Side Rail Extra Height	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ Patent Angle Indicator For Back-rest Position	2pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Trapeze Socket
☐ Oxygen Cylinder Holder
☐ Basket

V8w
Electric Bed
SKD-C

75°±10°
Back-rest

46°±10°
Knee-rest

290mm
Height
Adjustment
455-745mm

12°±2°
Reverse
Trendelenburg

12°±2°
Trendelenburg

LINAK
Motor



Technical Parameters

External Length	2085mm
External Width	1020mm
Bed Height Min.	455mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg


Technical Configuration

✓ LINAK Motor	4pcs
✓ Hand Controller	1pc
✓ 5" Hardshell Castors	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Trapeze Socket
☐ Basket

V4w
Manual Bed
SK-CF1



Technical Parameters

External Length	2130mm
External Width	1020mm
Bed Height Min.	475mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Cranks	4pcs
✓ 5" Hardshell Castors	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Trapeze Socket
☐ Basket

V3v
Manual Bed
SK-CD1



Technical Parameters

External Length	2130mm
External Width	1025mm
Bed Height Min.	430mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Cranks	3pcs
✓ 5" Hardshell Castors	4pcs
✓ PP Side Rail Extra Height	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ Patent Angle Indicator For Back-rest Position	2pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Utility Shelf
☐ Trapeze Bracket/Pole

V2k
Manual Bed
SK-C1



Technical Parameters

External Length	2130mm
External Width	995mm
Bed Height	500mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

☒ Manual Crank	2pcs	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Monkey Pole	
☒ PP Bed Ends with Safe Lock with Φ 65 Bumper Wheels	1 set	☐ Trapeze Socket	
☒ PP Bed Platform with Integrated Mattress Retainer	1 set	☐ Basket	
☒ IV Pole 4 Hooks	1 set		

V6k
Electric Bed
SKD-D



Technical Parameters

External Length	2085mm
External Width	995mm
Bed Height Min.	420mm
Bed Height Max.	730mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

☒ Linear Motor	3pcs	☒ IV Pole 4 Hooks	1 set
☒ Hand Controller	1 pc	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Monkey Pole	
☒ PP Bed Ends with Safe Lock with Φ 65 Bumper Wheels	1 set	☐ Trapeze Socket	
☒ PP Bed Platform with Integrated Mattress Retainer	1 set	☐ Basket	



V3k
Manual Bed
SK-CD1



Technical Parameters

External Length	2130mm
External Width	995mm
Bed Height Min.	420mm
Bed Height Max.	730mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

☒ Manual Crank	3pcs	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Monkey Pole	
☒ PP Bed Ends with Safe Lock with Φ 65 Bumper Wheels	1 set	☐ Trapeze Socket	
☒ PP Bed Platform with Integrated Mattress Retainer	1 set	☐ Basket	
☒ IV Pole 4 Hooks	1 set		



K6k
Electric Bed
SKD-D

84°±10°
Back-rest

47°±10°
Knee-rest

310mm
Height
Adjustment
400-710mm

Technical Parameters

External Length	2105mm
External Width	980mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1925×820mm
Safe Working Load	200kg

Technical Configuration

☒ Linear Motor	3pcs	☒ IV Pole Prevision	4pcs
☒ Hand Controller	1 pc	☒ Drainage Hook	4pcs
☒ 5" Covered Castors	4pcs	☒ PP Mattress Retainer	2pcs
☒ Foldable Side Rail	1 set	☐ Basket	
☒ PP Bed Ends with Safe Lock	1 set		
☒ IV Pole 4 Hooks	1 set		

K3k
Manual Bed
SK-CD1

85°±10°
Back-rest

45°±10°
Knee-rest

310mm
Height
Adjustment
400-710mm

Technical Parameters

External Length	2160mm
External Width	980mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1930×830mm
Safe Working Load	200kg

Technical Configuration

☒ Manual Crank	3pcs	☒ Drainage Hook	2pcs
☒ 5" Covered Castors	4pcs	☐ Basket	
☒ Foldable Side Rail	1 set		
☒ PP Bed Ends with Safe Lock	1 set		
☒ IV Pole 4 Hooks	1 set		
☒ IV Pole Prevision	4pcs		

K2k
Manual Bed
SK-C1

80°±10°
Back-rest

35°±10°
Knee-rest

Technical Parameters

External Length	2160mm
External Width	980mm
Bed Height	500mm
Mattress Platform	1930×830mm
Safe Working Load	200kg


Technical Configuration

☒ Manual Crank	2pcs	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Basket	
☒ PP Bed Ends with Safe Lock	1 set		
☒ IV Pole 4 Hooks	1 set		



SK-HD03

Rehabilitation Bed



Maximum patient weight

The extra stable and sturdy frame allows a safe load of up to 320 kg.



SK-HD01

Electric Homecare Bed



SK-HD02

Electric Homecare Bed



SK-HD04

Electric Homecare Bed




SK-HD05

Electric Homecare Bed



Model	SK-HD01	SK-HD02	SK-HD03	SK-HD04	SK-HD05
External size	2140×1055mm	2120×1025mm	2195×1345mm	2060×1050mm	2060×1070mm
Mattress platform	2000×900mm	2000×900mm	2010×1200mm	2000×900mm	2000×900mm
Back-rest adjustment	0-72°	0-72°	0-70°	0-68°	0-68°
Leg-rest adjustment	0-28°	0-28°	0-24°	0-28°	0-28°
Reverse trendelenburg	0-12°	0-12°	0-12°	0-15°	0-17°
Trendelenburg	0-12°	0-12°	0-12°	0-18°	0-15°
Lowest & highest position	190-610mm	390-810mm	230-660mm	260-760mm	300-800mm
Safe working load	215kg	215kg	320kg	215kg	215kg
Mattress size	L2000×900mm/800mm	L2000×900mm/800mm	L2000×1200mm	L2000×900mm/800mm	L2000×900mm/800mm
Guardrail	Wooden full-length	Wooden full-length	Wooden full-length	Matel collapsible	Wooden split
Construction	Foldable	Foldable	Foldable	Non-foldable	Non-foldable
Feature	Super low position	Economy model	Super load bearing Super wide width	Enormous height range	Enormous height range

Q5k
Electric Children Bed
SKD-A



360° View


Technical Parameters

External Length	1930mm
External Width	880mm
External Height	600mm
Mattress Platform	1710×780mm

Technical Configuration

✓ Linear Motor	2pcs	✓ 5" Double Side Castors	4pcs
✓ Four-fold Aluminum Alloy Side Rails	1set	✓ IV Pole Prevision	4pcs
✓ Bed End With Medical Record Card Holder	1set	✓ IV Pole	1set
✓ Detachable PP Bed Surface With Mattress Holder	1set	✓ Drainage Hook	2pcs
✓ Hand Controller	1pc		


Q2w
Children Bed
SK-C1



Technical Parameters

External Length	2185mm
External Width	880mm
External Height	540mm
Mattress Platform	1940×780mm

CQ2k
Children Bed
SK-C2



Technical Parameters

External Length	1975mm
External Width	880mm
External Height	540mm
Mattress Platform	1730×780mm



Technical Configuration


✓ Manual Cranks	2pcs	✓ IV Pole Prevision	4pcs
✓ Six-fold Aluminum Alloy Side	1set	✓ IV Pole	1set
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Drainage Hook	2pcs
✓ Detachable PP Bed Surface With Mattress Holder	1set	✓ Shoes Holder	1set
✓ Bed End With Medical Record Card Holder	1set		



Technical Configuration

✓ Manual Cranks	2pcs	✓ 5" Covered Castors	1set
✓ Bed End With Medical Record Card Holder	4pcs	✓ IV Pole Prevision	4pcs
✓ Twelve-fold Foldable Stainless Steel Aluminum Alloy Side Rails	1set	✓ IV Pole	1set
✓ Detachable PP Bed Board With Built-in Mattress Holder	1set	✓ Drainage Hook	2pcs

CQ8k
Electric Children Bed
SKD-C



360° View


Technical Parameters

External Length	1975mm
External Width	880mm
Bed Height Min.	480mm
Bed Height Max.	770mm
Mattress Platform	1730×780mm

Technical Configuration

✓ LINAK Motor	4pcs	✓ PP Bed Ends with Safe Lock	1 set
✓ Hand Controller	1 pc	✓ IV Pole	1 set
✓ PP Platform with Integrated Mattress Retainer	1 set	✓ IV Pole Prevision	4pcs
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Side Rail 12 Columns with Extra Protection	1 set	✓ Electric CPR	

SK-CHD01
Super Low Children Bed




Technical Parameters

External Length	2040mm
External Width	926mm
Bed Height Min.	268mm
Bed Height Max.	688mm
Mattress Platform	1900×780mm

Technical Configuration

✓ Linear Motor	4pcs
✓ Wooden Side Rails	1 set
✓ Head & Foot Board	1 set
✓ Castors	4pcs
✓ Remote Controller	1 pc

CD3q
Children Bed
SK-CD2



Technical Parameters

External Length	2030mm
External Width	1045mm
Bed Height Min.	500mm
Bed Height Max.	770mm
Mattress Platform	1705×900mm

Technical Configuration

✓ Manual Cranks	3pcs	✓ IV Pole	1 pc
✓ 5" Hardshell Castors	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Full Length Side Rail	1 set		
✓ ABS Grace Bed Ends with Safe Lock	1 set		
✓ IV Pole Prevision	4pcs		

CX2x
Children Bed
SK-C2



Technical Parameters

External Length	1910mm
External Width	930mm
Bed Height	580mm
External Height	1260mm
Mattress Platform	1775×710mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole	1pc
✓ 5" Double Side Castors	4pcs	✓ Buffer Air Spring	4pcs
✓ Full Length Side Rail	1set		
✓ Plastic Bed Ends	1set		
✓ IV Pole Prevision	4pcs		

CR0q
Children Bed
SK-A2



Technical Parameters

External Length	1600mm
External Width	750mm
Bed Height	500mm
External Height	1040mm
Mattress Platform	1580×660mm

Technical Configuration

✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Full Length Side Rail	1set		
✓ Integrated Bed Ends	1set		
✓ IV Pole	1pc		
✓ IV Pole Prevision	4pcs		

CR2q
Children Bed
SK-C2



Technical Parameters

External Length	1800mm
External Width	860mm
Bed Height	600mm
External Height	1110mm
Mattress Platform	1600×680mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole Prevision	4pcs
✓ Full Length Side Rail	1set	✓ Drainage Hook	2pcs
✓ Intergrated Bed Ends	1set	✓ Shoes Holder	1pc
✓ Anti-slip Bed Feet Cover	4pcs		
✓ IV Pole	1pc		

X01-1

Baby Crib
SK-AF3



Technical Parameters

External Length	880mm
External Width	520mm
Bed Height Min.	770mm
Bed Height Max.	940mm

Technical Configuration

✓ Gas Spring	2pcs
✓ Bassinet	1pc
✓ Φ 100 Double Side Castors	4pcs
✓ 2cm Thickness Mattress	1pc

X01-5

Baby Crib
SK-AE3



Technical Parameters

External Length	810mm
External Width	520mm
Bed Height	990mm

Technical Configuration

✓ Bassinet	1pc	✓ IV Pole	1set	✓ 2cm Thickness Mattress	1pc
✓ 3" Covered Castors	4pcs	✓ Utility Basket	1pc		
✓ IV Pole Prevision	2pcs	✓ Adjustable Air Spring	1pc		

X01

Baby Crib
SK-AE3



Technical Parameters

External Length	810mm
External Width	490mm
Bed Height	960mm

Technical Configuration

✓ Bassinet	1pc	✓ IV Pole	1set
✓ 3" Castor With Brake	4pcs	✓ 2cm Thickness Mattress	1pc
✓ IV Pole Prevision	1pc	✓ Stainless Steel Frame	

X20

Pediatric Exam Table
SKZ-A

Technical Parameters

External Length	1800mm
External Width	650mm
Bed Height	830mm

Technical Configuration

✓ Material:laminate Wood Structure With PU Mattress
✓ Foot Part Can Be Slided Out For Tall Childred
✓ With Lockers At Table Side
✓ With Cartoon Stickers



SK-PB025

3D Swing Spinal Therapy Bed



Technical Parameters

External Length	1700mm
External Width	580mm
Bed Height Min.	570mm
Bed Height Max.	890mm
Safe Working Load	200kg

Technical Configuration

Adjustable width headrests	2pcs	Tension controlled drop devices	4pcs	High-resistance springs	2pcs
Fixed armrests	2pcs	Stainless steel handles	4pcs	Stainless steel operating lever	1pc
Tissue holder	1set	Bakelite adjustment knob	1set	Leather bed surface	1set
Electric actuator	3pcs	Large hydraulic gas rod	1pc	Square tube bed frame	1set
Foot controller	1set	Small hydraulic gas rods	2pcs	Foot restraint strap	1set

SK-PB023

Electric Spinal Therapy Bed



Technical Parameters

External Length	1800mm
External Width	550mm
Bed Height Min.	520mm
Bed Height Max.	740mm
Safe Working Load	200kg

Technical Configuration

Adjustable width headrests	2pcs	Stainless steel handles	4pcs
Fixed armrests	2pcs	Bakelite adjustment knob	1set
Tissue holder	1pc	Leather bed surface	1set
Electric actuator	3pcs	Square tube bed frame	1set
Foot controller	1set	3-inch castors	2pcs
Tension controlled drop devices	4pcs	Anti-slip foot pads	2pcs

SK-PB021

Electric Spinal Therapy Bed



Technical Parameters

External Length	1780mm
External Width	550mm
Bed Height Min.	550mm
Bed Height Max.	780mm
Safe Working Load	200kg

Technical Configuration

Adjustable width headrests	2pcs	Square pillow	1pc
Fixed armrests	2pcs	Electric actuator	2pcs
Hydraulic gas rod	2pcs	Foot controller	1set
Tension controlled drop devices	5pcs	PVC leather bed surface	1set
Foot pedals	8pcs	Square tube bed frame	1set
Bakelite adjustment knob	1set	Anti-slip foot pads	4pcs

SK-PB018

Spinal Therapy Bed



Technical Parameters

External Length	1800mm
External Width	550mm
Bed Height	530mm
Safe Working Load	200kg

Technical Configuration

Adjustable width headrests	2pcs	Bakelite adjustment knob	1set
Fixed armrests	2pcs	Durable PVC leather bed surface	1set
Hydraulic gas rod	1pc	Low-carbon square tube bed frame	1set
Tension controlled drop devices	3pcs	3-inch wear-resistant nylon casters	2pcs
Foot pedals	6pcs	Anti-slip foot pads	4pcs

X36-1
Electric Examination Table




Technical Parameters

External Length	1960mm
External Width	680mm
Table Height Min.	470mm
Table Height Max.	910mm

Technical Configuration

✓ Adjustable arm support	1pc	✓ Brake Pedal	4pcs
✓ Motor	1pc	✓ Hand Controller	1pc
✓ Gas spring	2pcs	✓ Manual Operator	1pc
✓ 3" Castor	4pcs		

X36
Electric Examination Table




Technical Parameters

External Length	1890mm
External Width	680mm
Table Height Min.	470mm
Table Height Max.	910mm

Technical Configuration

✓ Pillow	1pc	✓ Brake Pedal	4pcs
✓ Motor	1pc	✓ Hand Controller	1pc
✓ Gas spring	2pcs		
✓ 3" Castor	4pcs		

X35
Electric Examination Table




Technical Parameters

External Length	1845mm
External Width	680mm
Table Height Min.	550mm
Table Height Max.	860mm

Technical Configuration

✓ Brake Pedal	2pcs	✓ Face Pillow	1pc
✓ Motor	1pc	✓ Tissue Holder	1pc
✓ Hand Controller	1pc		
✓ 2" Casters	4pcs		

X14
Hydraulic Examination Table
SKZ-BD



Technical Parameters

External Length	1940mm
External Width	700mm
Table Height Min.	490mm
Table Height Max.	890mm

Technical Configuration

✓ Hydraulic Pump	1pc	✓ Toilet Paper Holder	1pc
✓ Adjustable Air Spring	2pcs		
✓ 5" Covered Castor	4pcs		
✓ Hydraulic pedal	2sets		

X07a
Examination Table
SKZ-A



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

☒ Anti-slip Bed Feet Cover	4pcs

X09
Examination Table
SKZ-B




Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

☒ Anti-slip Bed Feet Cover	4pcs
☒ Mechanical Lever to Adjust Back Rest	1set

A048
Gynecological Examination Table
SKZ-B




Technical Parameters

External Length	1810mm
External Width	690mm
Table Height	880mm

Technical Configuration

☒ Drawer	3pcs	☒ Lithotomy Rod	1set
☒ Step Stool	1pc	☒ Stainless Steel Basin	1pc
☒ Retractable Foot Board	1pc	☒ Paper Roll Holder	1pc

A98-8
Electric Obstetric Bed




Technical Parameters

External Length	2300mm
External Width	1050mm
Table Height Min.	480mm
Table Height Max.	880mm

Technical Configuration

✓ Guardrail	2pcs	✓ Delivery Holder	1pc	✓ Mattress	1set
✓ Head&foot Board	1pc	✓ Foot Support	1pc	✓ Basin	1pc
✓ Motor	1set	✓ 6" Castor, Central Locking System	4pcs	✓ IV Pole	1pc
✓ Obstetric Handle	2pcs	✓ Wall Bumper Wheel	2pcs	✓ Backup Battery	1pc
✓ Auxiliary Table	1pc	✓ Night Light	1pc		

A98-3
Electric Obstetric Bed



Technical Parameters

External Length	2120mm
External Width	1120mm
Table Height Min.	595mm
Table Height Max.	905mm

Technical Configuration

✓ Pull Handle	2pcs	✓ Stainless Steel Basin	1pc
✓ Leg support	2pcs	✓ Instrution Manual	1pc
✓ Foot support	2pcs	✓ 6" Castor	1pc
✓ Hand Remote Controller	1pc	✓ Mattress	1set

A98-4

Electric Obstetric Bed

360mm
Height
Adjustment

70°
Back-rest

10°±2°
Reverse
Trendelenburg

10°±2°
Trendelenburg

65°/5°
Leg Plate
up/down



Technical Parameters

External Length	2050mm
External Width	960mm
Table Height Min.	640mm
Table Height Max.	1000mm

Technical Configuration

✓ Linear Motor	3pcs	✓ Central Locking Pedal	2pcs	✓ Hand Controller	1pc
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Hand Grip	1 pair	✓ Power Line	1pc
✓ ABS Side Rail With Angle Indicator	1 set	✓ Leg Holder	1 pair	✓ Leg Strip	1 Pair
✓ ABS Head Board With Safe Lock	1pc	✓ Basin	1pc	✓ Iv Pole Prevision	2pcs

A99-13

Electric Obstetric Bed

250mm
Height
Adjustment

50°
Back Plate
up

22°
Reverse
Trendelenburg

22°
Trendelenburg

90°
Leg Plate
Swing Out

16°
Seat Plate
up



Technical Parameters

External Length	1860mm
External Width	600mm
Table Height Min.	730mm
Table Height Max.	980mm

Technical Configuration

✓ Armrest	2pcs	✓ SS Rinsing Basin	1pc
✓ Knee Curtch	2pcs	✓ Pedal Controller	1pc
✓ Pull Handle	2pcs	✓ LED Examination Light	1pc
✓ Anesthesia Frame	1pc		

A99-7

Electric Obstetric Bed

250mm
Height
Adjustment

55°/18°
Back Plate
up/down

35°/3°
Seat Plate
up/down

90°
Leg Plate
Swing Out



Technical Parameters

External Length	1980mm
External Width	720mm
Table Height Min.	630mm
Table Height Max.	880mm

Technical Configuration

✓ Bed Frame	1pc	✓ Hand Bars	2pcs	✓ Castor	1pc
✓ Foot Pedal	1pc	✓ Leg Support with Cushion	2pcs		
✓ Hand Controller	1pc	✓ Arm Support with Cushion	2pcs		
✓ Anesthesia Screen	1pc	✓ Stainless Steel Basin	1pc		

A99-5
Electric Obstetric Bed

300mm
Height
Adjustment

55°/10°
Back Plate
up/down

20°/10°
Seat Plate
up/down

90°
Leg Plate
Swing Out



Technical Parameters

External Length	1985mm
External Width	665mm
Table Height Min.	675mm
Table Height Max.	975mm

Technical Configuration

✓ Arm Support	2pcs	✓ Tightener	2pcs	✓ Slider	4pcs
✓ Handle	2pcs	✓ Hand Controller	1pc		
✓ Leg Support	2pcs	✓ Foot Switch	1pc		
✓ Stirrup	2pcs	✓ Stainless Steel Basin	1pc		

A105
Manual Obstetric Table

250mm
Height
Adjustment

50°/10°
Back Plate
up/down

20°
Reverse
Trendelenburg

10°
Trendelenburg

90°
Leg Plate
Swing Out



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	700mm
Table Height Max.	950mm

Technical Configuration

✓ Bed Frame	1pc	✓ Leg Support with Cushion	2pcs
✓ Mattress	1pc	✓ Arm Support with Cushion	2pcs
✓ Anesthesia Screen Frame	1pc	✓ Stainless Steel Basin	1pc
✓ Hand Bars	2pcs		

A99-8
Electric Gynecological Exam Couch

70°/20°
Back Plate
up/down

40°
Seat Plate
up

490mm
Height
Adjustment



Technical Parameters

External Length	1400mm
External Width	740mm
Table Height Min.	610mm
Table Height Max.	1100mm

Technical Configuration

✓ Foot Switch	1pc	✓ Removable Handles	2pcs
✓ Power Switch	1pc	✓ Pull-out Stainless Steel Basin	1pc
✓ Seamless Foaming Mattress	1set	✓ 5" Double Sided Castor	1set
✓ Leg Support	2pcs		

A045-1
Gynecological Table
SKZ-C

70°±10°
Back Plate
up

90°±10°
Leg Plate
down



Technical Parameters

External Length	1820mm
External Width	590mm
Table Height	800mm

Technical Configuration

✓ Leg Support	2pcs	✓ Mechanical Lever	2sets
✓ Hand Bar	1set	✓ Stainless Steel Basin	1pc
✓ Anti-slip Bed Feet Cover	4pcs		



A106-5
Electric-Hydraulic
Operating Table

The A106-5 operating table is an electro-hydraulic operating table with a modular table top and unified interfaces. The head board, leg board and back board can be assembled and matched arbitrarily according to the surgical requirements. The head and leg boards are interchangeable. After the exchange, the system can be Automatically identify the direction of the head and leg boards.

30°
Reverse
Trendelenburg

30°
Trendelenburg

22°
Lateral Tilt
left/right

50°/90°
Head Plate
up/down

85°/45°
Back Plate
up/down

80°/90°
Leg Plate
up/down

100mm
Kidney
Bridge Up

300mm
Horizontal
Sliding
C type Arm

110°
Flex

220°
Reflex



Operating Table

Technical Parameters

External Length	2150mm
External Width	555mm
Table Height Min.	600mm
Table Height Max.	1100mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair		
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		



A302
Electric Operating Table
SKL-A

A302 is a highly automated electric Multi-functional operating table. Its elevation, lowering, trendelenburg, reverse trendelenburg, lateral tilt, back folding, vertical translation and leg folding are all electrically operated. The operating table is equipped with a digital display hand controller to make operation easy and intuitive. With Imported large diameter castors, smooth movement, stable brake system. The standard dual control system guarantees the safety of operation.

24°
Reverse
Trendelenburg

24°
Trendelenburg

23°
Lateral Tilt
left/right

33°/90°
Head Plate
up/down

68°/22°
Back Plate
up/down

18°/85°
Leg Plate
up/down

110mm
Kidney
Bridge Up

350mm
Horizontal
Sliding
C type Arm

100°
Flex

220°
Reflex



Technical Parameters

External Length	2000mm
External Width	555mm
Table Height Min.	690mm
Table Height Max.	990mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		



A301
Electric Operating Table
SKL-A

A301 is a multi-functional electric operating table. Its mattress is made of polyurethane molded material, can be sterilized and cleaned at high temperature, and there are no gaps so to avoid bacterial growth. The platform is made of high-purity phenolic resin board applicable for C-arm, and can be slid up to 350mm in one direction. The X-ray photography can carry out the whole body radiography examination without moving the patient, no blind spots or shadows.

24°
Reverse
Trendelenburg

24°
Trendelenburg

23°
Lateral Tilt
left/right

30°/90°
Head Plate
up/down

66°/22°
Back Plate
up/down

90°
Leg Plate
down

90°
Leg Plate
Horizontal
Positioning

110mm
Kidney
Bridge Up

350mm
Horizontal
Sliding
C type Arm



Operating Table

Technical Parameters

External Length	2000mm
External Width	555mm
Table Height Min.	730mm
Table Height Max.	980mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		

A100-4A

Electric Operating Table

SKL-A

A100-4A is an electric operating table with super load-bearing and extra-long sliding distance , which can be applicable for C-arm. The built-in large-capacity battery can continue to support all operations in the event of a power failure. And the operating table adopts a dual-control system to ensure the safety of the operation.



Technical Parameters

External Length	2115mm
External Width	550mm
Table Height Min.	710mm
Table Height Max.	1010mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports (with straps)	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests (with straps)	1 pair		
☒ Lumbar Supports	1 pair		

A100-4

Electric Operation Table

SKL-A

A100-4 is an electric operating table that can meet the needs of various surgeries. It is driven by low-voltage DC motor, with safe and stable performance. Super wide table top, super long sliding distance. It can be applicable for C-arm. It has a high level of automation, low noise, strong reliability, and imported original parts for key parts. It is a perfect and ideal operating table.



Technical Parameters

External Length	2080mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm


Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports (with straps)	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests (with straps)	1 pair		
☒ Lumbar Supports	1 pair		

A501

Electric-Hydraulic Operating Table

A501 is an advanced, multi-functional electric-hydraulic operating table. Its table top is made of high-strength synthetic board and can transmit X-rays. In addition, the platform can be equipped with a 180° rotation function to realize head and leg interchange, which is applicable for C-arm perspective and placement of some special positions.



Technical Parameters

External Length	2125mm
External Width	550mm
Table Height Min.	630mm
Table Height Max.	980mm


Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair		
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		

A2000

Electric Operating Table SKL-A

A2000 is a basic electric operating table with a very high cost effective, it can fulfill the requirement of general surgeries. Its built-in battery can continue to support all operations in the event of a power failure.




Technical Parameters

External Length	2050mm
External Width	480mm
Table Height Min.	680mm
Table Height Max.	930mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		

AM052C-1
Manual Operating Table
SK1F-R



Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min.	790mm
Table Height Max.	1040mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports (with straps)	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests (with straps)	1 pair		
☒ Lumbar Supports	1 pair		

Multi-functional Orthopedic Frame

Gel Mattress

Carbon Fiber Multi-functional Orthopedic Frame

Upper Limb Traction Device

Multi-purpose Head Rack

Shoulder Traction Frame

Stirrup Pneumatic Power Leg Rack

Bow Spine Operating Supporter

Mechanical Power IV Pole

Operating Table Accessories Cart

Knee Surgery Rack

Double Table Top

Prone Body Supporting Frame

Upper Limb Operating Table

Bed Surface Widening Module

External Kidney Bridge

Dirt Basin

Body Safe Strap



SHADOW MANAGEMENT

Everything has been designed to enable the medical team to focus exclusively on the surgical operation thanks to the perfect shadow dilution. A large number of individual divergent light beams create a complete light patch, ensuring that the light field will remain circular and shadow-free. Alno Light20 provides the best illumination, even when surgeons are blocking the LED.



Bulb

OSRAM LED bulbs with 72 beads and 24 triangular lenses per head, providing long-lasting, bright, and uniform illumination.



Color Temperature&illumination

Advanced LED cold light with up to 160,000 Lux, 10-level adjustment, and 3500-5500K color temperature for soft, non-glare illumination.



5-inch Extra-large LCD

5-inch large LCD with simple interface, 30° adjustable angle, and seamless design for easy disinfection.

	Hours	Kelvin	Angle Adjustment
OSRAM	60000	3500-5500	±30°

The head of the lamp is windmill-shaped, Saikang's unique (appearance patent) design. The back is equipped with a soft ambient light. The top offers a soft ambient light, simple and stylish while also being useful for low-light surgery (laparoscopic). Exceeding the infection control standards of modern operating rooms.



User friendly for multiple surgeons and nurses at the same time



Head insert design: No heating over the surgeons head (max. Allow is 1C°)

Alno /series
LED Operation Lamp



Operating Handle

Detachable PPSU operating handle, heat- and pressure-resistant for safe steam sterilization.



Intelligent Shadow Management

Smart sensors adjust lighting to eliminate shadows and maintain clear visibility.

Lighting Mode

It has three lighting modes (standard, endoscopic, and low-light), as well as red light compensation, green light compensation, and shadow compensation.



STANDARD



ENDOSCOPIC



LOW-LIGHT

Alno Light20



Alno Light10



Alno Light M



	Alno Light20 SK-SD700/700	Alno Light10 SK-SD700	Alno Light M
Central illuminance (Ec)	40000-160000 lx (10 levels)	40000-160000 lx (10 levels)	40000-160000 lx (10 levels)
Spot diameter (d10)	200-300mm (5 levels)	200-300mm (5 levels)	200-300mm (5 levels)
Spot distribution (d50/d10)	55%	55%	55%
Depth of the light column (60% of Ec)	700mm	700mm	700mm
Color temperature (Tc)	3500-5500K (11 levels)	3500-5500K (11 levels)	3500-5500K (11 levels)
Colour rendering index (Ra)	98	98	98
Maximum total irradiance (Ee)	≤550 W/m²	≤550 W/m²	≤550 W/m²
Radiant energy (Ee/Ec)	≤3.5mW/(m².lx)	≤3.5mW/(m².lx)	≤3.5mW/(m².lx)



Alnno X71

Medical Bridge Ceiling Pendant
(separate wet and dry areas)

Modern ceiling pendant with separate wet & dry areas, offering flexible platforms, drawers, and infusion holders for efficient and safe medical rescue.



Zoned Lighting System

Smart zoned ICU lighting with adjustable brightness and color; all controlled by touch panel.



Information Display

Built-in bed info display for quick patient data access and improved efficiency.

Alno / series

Suitable for ICU, NICU, PICU, EICU room

Wet and Dry Areas Separated

Separated wet and dry areas ensure clear layout, minimize interference, and enhance safety.



Translation System

Translation system with 500mm range, smooth control, and precise positioning for multi-directional care.



Suspension System

Both the tower and arms can rotate 340°, with built-in limit devices to prevent drift.

Rotate
340°



Telescopic Extension Arms

Telescopic extension arms with high load capacity, flexibly adapt to infusion needs and optimize space.



Braking System

One-button remote braking with indicator lights ensures quick, precise control.



Drawer

Durable 30kg drawer with smooth silent rails and optional sensor light for convenience.



Cable Management



Cable Hub



Stainless steel basket (Optional)



Suction tube box (Optional)



Drawer divider box (Optional)



Computer keyboard drawer (Optional)



Flashlight + flashlight placement box (Optional)



Stethoscope + stethoscope hanger (Optional)

Alnno L21

Electric Double-Arm Pendant

- Lightweight, enclosed design – safe and flexible.
- Electric double-arm – wide treatment coverage.
- Supports multiple devices – space-efficient.
- Boosts OR efficiency – saves space.



Alnno X22

Double-Arm Pendant

- Lightweight, compact design – easy and safe to use.
- Double-arm movement – extends treatment reach.
- Accommodates multiple devices – maximizes workspace.
- Enhances OR workflow – saves space and time.




Mattress (OEM Size)




SKP028

Pressure Reliving Mattresses with Gel




SKP040

Cube Design Pressure-relieving Mattress



SKP041

Pressure-relieving Mattress





Mattress (OEM Size)

SKP011



SKP001



SKS009-2
Phenolic Resin Bedside Table

Technical Parameters

External Size	590x480x840mm
Dining Table Max. Height	1000mm

Technical Configuration

✓ Top Guardrail	Aluminum alloy
✓ Drawer	1 pc
✓ Cabinet	1 pc
✓ Casters	4 pcs
✓ Side Pull-up Writing Board	1 pc

SKS02-W
ABS Bedside Table
475×465×795mm(L×W×H)



SKS03-W
ABS Bedside Table
480×450×950mm(L×W×H)



SKS007-4
Steel Bedside Table
440×460×810mm(L×W×H)



SKS036
Bedside Table
440×440×790mm(L×W×H)



SKS008
ABS+Steel Bedside Table
500×510×705mm(L×W×H)



SKS013-1
ABS+Steel Bedside Table
575×500×750mm(L×W×H)





SKS025
Phenolic Resin Bedside Table
475×475×800mm(L×W×H)

SKH046-11
Overbed Table

- Technical Parameters
800×395×860-1160mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ Phenolic Resin Table Top



SKH046-2
Overbed Table



- Technical Parameters
890×445×760-1050mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ PE Table Top with Cup Socket



SKH046-14
Overbed Table



- Technical Parameters
890×400×855-1155mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ ABS Table Top with Cup Socket and Push Handle



SKH042
Overbed Table



- Technical Parameters
800×410×810-1110mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ MDF Table Top



SKH046-SC
Telescopic Dining Board



- Technical Parameters
890-1130×350mm
- Technical Configuration
- ✓ Hook



SKE-136
Electric Dalysis Chair
SK1D-0F2



Technical Parameters

External Length	2020mm
External Width	600mm
Bed Height Min.	620mm
Bed Height Max.	870mm

Technical Configuration

✓ Hand Controller	1pc	✓ Arm Support	2pcs
✓ USB Port	1pc	✓ Plug-in Headrest (with a Small Pillow)	1pc
✓ 3" Brake Wheel	4pcs	✓ Seamless Soft - Padded Seat Cushion	2pcs
✓ Small Dining Table	1pc		

SKE001

Attendant Chair

Technical Parameters

Seat Size (L×W×H)	680×720×940mm
Bed Size (L×W×H)	680×1900×590mm



SKE001-5

Therapy Chair

Technical Parameters

Seat Size (L×W×H)	780×650×1060mm
Bed Size (L×W×H)	2060×650×450mm



SKE001-16

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	1000×800×860mm
Bed Size (L×W×H)	1920×600×440mm



SKE001-19

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	1600×700×655mm
Bed Size (L×W×H)	2200×700×655mm



SKE091

Blood Donation Chair

Technical Parameters

Seat Size (L×W)	1750×600/940mm
Seat Height	480-540mm



SKE941

Therapy Chair

Technical Parameters

Seat Size (L×W×H)	940×770×1250mm
Bed Size (L×W×H)	1930×770×530mm



SKE001-10

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	850×1000×850mm
Bed Size (L×W×H)	600×1850×430mm



SKE014-9

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	1250×1000×850mm
Bed Size (L×W×H)	1030×1850×430mm





SKE909

Reclining Sofa

Technical Parameters

Seat Size (LxW)	640×800mm
Bed Size (LxW)	640×1950mm



SKE087

Reclining Chair

Technical Parameters

Back size (LxW)	750x540mm
Seat Size (LxW)	540x540mm
Leg Size (LxW)	560x540mm



SKE090

Blood Donation Chair

Technical Parameters

Seat Size (LxWxH)	910×480×960mm
-------------------	---------------


Technical Configuration

- ✓ Powder coated steel structure
- ✓ Upholstery PU with foam
- ✓ Padded arm rest, adjustable 360 degree swivel
- ✓ With cabinet at side for storage
- ✓ Rubber feet cushion




SKE011-1
Treat-waiting Chair

Material: PU
1630×650×865mm(LxWxH)




SKE011-6
Treat-waiting Chair

Material: PU
1655×640×820mm(LxWxH)



SKE011-8
Treat-waiting Chair

Material: PU
1780×600×820mm(LxWxH)





360° View

SK-CL001 SK-CL501
Curtain Track



SKH049-7
Phenolic Resin Screen
Full Size: 1780×1650mm
Folded Size: 594×412×175mm



SKH049
Stainless Steel Screen
Size: 2000×1870mm



DM-Z104
Office Table
Size: 1400×600-1040×750mm



DM-Y106
Office Chair
Size: 650×640×965mm



SKE801
Nurse Stool
Size: Φ360×450-570mm



SKE020-1
Step Stool
Size: 490×410×170/340mm



SKH034
Trash Can



12L

SKH034-5
Kick Bucket
Size: Φ400×H360×14L



14L

SKH041-3

I.V. Stand

Function

Height Adjustment 1350-2250mm

Material

IV Pole Outer Tube: Φ25×1000mm (Stainless Steel)

IV Pole Inner Tube: Φ16×1200mm (Stainless Steel)



SKH041(1)

I.V. Stand

Function

Height Adjustment: 1350-2250mm

Material

IV Pole Outer Tube

Φ25×1000mm (Stainless Steel)

IV Pole Inner Tube

Φ16×1200mm (Stainless Steel)



SKH041(18)

I.V. Stand

Function

Height Adjustment: 1200-2000mm

Material

IV Pole Outer Tube

Φ32 (Stainless Steel)

IV Pole Inner Tube

Φ19mm (Stainless Steel)



SKH041(16)

I.V. Stand

Function

Height Adjustment: 1500-2000mm

Material

IV Pole Outer Tube



Φ32mm (Stainless Steel)

IV Pole Inner Tube

Φ22-Φ25mm (Stainless Steel)



SKB1C10
Electric Stair Stretcher



Technical Parameters

Full Size	1225×510×1620mm
Folded Size	1080×510×270mm

SKB039(G)
Stretcher Trolley
SKC-3B10



Technical Parameters

Full Size	2000×560×1000mm	Knee-rest	20°
Folded Size	2000×560×380mm	Mattress thickness	PVC, 100 mm
Back-rest	75°	Load Bearing	200kg

SKB039(H)
Stretcher Trolley
SKC-3B4



Technical Parameters

Full Size	1950×570×910mm	Mattress thickness	PVC, 100 mm
Folded Size	2030×570×300mm	Load Bearing	180kg
Back-rest	75°		

SKB039(H)
Stretcher Trolley
SKC-3B4



Technical Parameters

Full Size	1950×570×910mm	Knee-rest	15°
Folded Size	2030×570×280mm	Mattress thickness	PVC, 100 mm
Back-rest	65°	Load Bearing	180kg

SKB037(A)
Stretcher Trolley
SKM-A
Size: 1960×600×810mm



SKB037(C)
Patient Trolley
SKM-G
Size: 1930×760×620-900mm



75°
Back-rest

35°
Knee-rest

280mm
Height
Adjustment
620-900mm

SKB038-4
Patient Trolley
SKM-E
Size: 1900×700×620-900mm



65°
Back-rest

280mm
Height
Adjustment
620-900mm

SKB037(B)
Patient Trolley
SKM-C
Size: 1960×520×770mm



70°
Back-rest

SKB038-2
Patient Trolley
SKM-E
Size: 1900×700×620-900mm



65°
Back-rest

280mm
Height
Adjustment
620-900mm

SKB041-3S
Hydraulic Patient Trolley
SKM-E
Size: 2150×790×635-895mm




60°
Back-rest

260mm
Height
Adjustment
635-895mm



SKB041-7
Patient Transportation Trolley
SKV-BE



Technical Parameters

Length	2200mm
Width	770mm
Height (X-ray Backrest)	550-860mm
Height (X-ray Platform)	600-910mm

Technical Parameters

☒ 80mm PU Mattress	1 pc	☒ Patient Restraints	2 pcs	☐ Foot-Board & Monitor Holder	1 pc
☒ 8" Castors, Central Locking	4 pcs	☒ PP Bottom Cover	1 pc	☐ Monitor Holder	1 pc
☒ Hydraulic Pumps	2 pcs	☒ 5-speed Foldable Guardrail	1 set	☐ Debris Basket	1 pc
☒ IV Pole	1 set	☒ Foot Pedals	1 set	☐ Weighting System	1 pc
☒ Back Angle Monitor	2 pcs	☒ 5th Directional Wheel	1 pc	☐ Pluggable Mobile IV pole	1 set
☒ Integrated Wall Bumper Wheels	4 pcs	☐ Oxygen Tanks Holder	1 pc		
☒ X-ray Translucent	1 set	☐ Paper Holder	1 pc		

SKB041-1

Patient Transportation Trolley
SKM-D



Technical Parameters

Length	1900mm
Width	640mm
Min. Height	550mm
Max. Height	850mm

Technical Configuration

☒ Crank System	1 set	☒ Directional Wheel	1 pc	☒ Central Locking Break Pedal	2 pcs
☒ Back Rest Gas Spring	1 pc	☒ IV Pole Previsions	2 pcs	☒ Transfer Mattress	1 pc
☒ ABS Side Rail	1 set	☒ IV Pole	1 set	☒ Directional Pedal	2 pcs
☒ 6" Castors with Central Locking	4 pcs	☒ Oxygen Cylinder Holder	1 pc		

SKB041-2

Patient Transportation Trolley
SKV-BC



Technical Parameters

Length	1910mm
Width	640mm
Min. Height	650mm
Max. Height	960mm

Technical Configuration

☒ Hydraulic Pump	2 pcs	☒ Directional wheel	1 pc	☒ Directional Pedal	2 pcs
☒ Back Rest Gas Spring	1 pc	☒ IV Pole Previsions	2 pcs	☒ Transfer Mattress	1 pc
☒ ABS Side Rail	1 set	☒ IV Pole	1 set		
☒ 8" Castors, Central Locking	4 pcs	☒ Central Locking Break Pedal	2 sets		

SKB041-6

Patient Transportation Trolley
SKM-F



Technical Parameters

Length	1930mm
Width	660mm
Min. Height	510mm
Max. Height	850mm

Technical Parameters

☒ Crank System	1 set	☒ Directional Wheel	1 pc	☒ Oxygen Cylinder Holder	1 pc
☒ Back Rest Air Spring	1 pc	☒ Directional Pedal	2 pcs	☒ Central Locking Break Pedal	4 pcs
☒ Plastic Side Rail	1 set	☒ IV Pole Previsions	4 pcs	☒ Transfer Mattress	1 pc
☒ 6" Castors, Central Locking	4 pcs	☒ IV Pole	1 set	☒ Safe Lock for Side Rail	2 pcs

SKB041-10

Patient Transportation Trolley
SKV-BE



Technical Parameters

Length	2150mm
Width	830mm
Min. Height	560mm
Max. Height	890mm

Technical Parameters

☒ Crank System	1 set	☒ Directional Wheel	1 pc	☒ Transfer Mattress	1 pc
☒ Back Rest Air Spring	2 pcs	☒ Directional Pedal	2 pcs	☒ Safe Lock for Side Rail	2 pcs
☒ Plastic Side Rail	1 set	☒ IV Pole Previsions	4 pcs	☒ Emergency Push Handle	1 set
☒ 8" Castors, Central Locking	4 pcs	☒ IV Pole	1 set	☐ Device platform	1 pc
☒ Hydraulic Pump	2 pcs	☒ Oxygen Cylinder Holder	1 pc		
☒ Wall Bumper Wheel	4 pcs	☒ Central Locking Break Pedal	4 pcs		

SKR-R10

ECG Machine Cart



SKR-R03

Laptop Cart



SKR-R04

Mobile Charting Cart



SKR-R07

Mobile Charting Cart



Model	SKR-R10	SKR-R03	SKR-R04	SKR-R07
Length	500mm (±10mm)	530mm (±10mm)	535mm (±10mm)	530mm (±10mm)
Width	500mm (±10mm)	520mm (±10mm)	520mm (±10mm)	620mm (±10mm)
Height	875mm (±10mm)	725-1045mm (±10mm)	725-1045mm (±10mm)	855-1100mm (±10mm)
Tabletop Size	335×365mm (±10mm)	450×445mm (±10mm)	450×505mm (±10mm)	480×460mm (±10mm)
Auxiliary Workbench	/	✓	/	/
Keyboard Case	/	✓	/	✓
Mouse Hanger	/	/	/	✓
Storage Box	/	/	/	✓
Column Material	Aluminum Alloy	Aluminum Alloy	Aluminum Alloy	Aluminum Alloy
Power System	/	Pneumatic Lifting	Pneumatic Lifting	Pneumatic Lifting
Base	Carbon Steel Base + ABS Cover	Cast Aluminum with Spraying Coating	Cast Aluminum with Spraying Coating	Cast Aluminum with Spraying Coating
Castor	5×2.5-inch double-sided castors, 2 of which are with brakes	4×4-inch double-sided castors, 2 of which are with brakes	4×4-inch double-sided castors, 2 of which are with brakes	4×4-inch double-sided castors, 2 of which are with brakes
Battery	/	/	/	Optional



SK-ET750

Emergency Trolley
Size: 750×530×1050mm

Technical Configuration

✓ Defibrillator shelf	1pc	✓ Castor	4pcs	✓ CPR Board	1pc
✓ Sliding Side Shelf	1pc	✓ Trash Can	2pcs	✓ Oxygen Bottle Holder	1pc
✓ Dust Baskets	1pc	✓ Central Lock	1pc		
✓ IV pole	2pc	✓ External Connectors	1 set		



SKR054-ET77A
Emergency Trolley
Size: 750×475×950mm



SKR054-ET
Emergency Trolley
Size: 760×475×945mm

SK-ET75077A
Emergency Trolley
Size: 750×475×950mm



SK-CT75077C2
Treatment Trolley
Size: 750×475×950mm



SK-AT75077A2
Anesthesia Trolley
Size: 750×475×950mm



SKR051-ET
Emergency Trolley
Size: 760×475×945mm



SKR054-AT
Anesthesia Trolley
Size: 760×475×945mm



SKR054-MT
Medicine Trolley
Size: 750×475×920mm



SKR058-ET

Emergency Trolley
Size: 790×600×1000/1100mm



SKR058-MT

Anesthesia Trolley
Size: 790×600×940/1020mm



SKR058-NT

Nursing Trolley
Size: 620×600×1000/1100mm



SKH040

Laundry Trolley
Size: 800×600×840mm



SKH040-1G

Laundry Trolley
Size: 800×600×865mm



SKH027

Laundry Trolley
Size: 900×450×1050mm



SKH006C

Instrument Trolley
(Disassembly structure)
Size: 670×450×900mm



SKH001

Treatment Trolley
Size: 670×450×900mm



SKH004

Instrument Trolley
Size: 670×450×900mm



SKH005

Fan-shaped Instrument Trolley
Size: L: 1550×460×900mm




SKH036


Stainless Steel Washing Sink
Size: 1200(W)×600(D)×1100(H)mm



Size: 600(W)×600(D)
×1100(H)mm



Size: 1800(W)×600(D)
×1100(H)mm



Size: 2400(W)×600(D)
×1100(H)mm



Technical Configuration

✓ Defibrillator shelf	1pc	✓ IV pole	1pc	✓ Sealed Handle	1pc
✓ Sliding Side Shelf	1pc	✓ Castor	4pcs		
✓ Writing Board	1pc	✓ Oxygen Bottle Holder	1pc		
✓ CPR Board	1pc	✓ Trash Can	1pc		

24hrs: 0086-18021231901 

www.saikangmedical.com

SAIKANG (Portugal Showroom)

Av. da Liberdade, nº 101 8150-101
São Brás de Alportel Portugal
TEL: +351 289 090 212

SAIKANG (Brazil Showroom)

Rua Lauro Muller 116, Sala 1503 - Botafogo,
Rio de Janeiro - RJ, 22290-160
E-mail: mehecobrasil@gmail.com

SAIKANG (China Headquarter)

No.35 Lehong Road, Zhangjiagang City,
Jiangsu Province, China.
TEL: +86 18021231901
E-mail: export@saikangmedical.com

SAIKANG (South East Asia Branch)

SAIKANG Medical (Thailand) Co., Ltd.
Ample Tower 909 Room 11/3 Debaratana Road,
Bangna Nuea, Bangkok 10260 Thailand
TEL: +66-936-569469

SAIKANG (US Branch)

Healthward Technology, INC
17800 castleton st ste 180, Los Angeles
TEL: +1 (603) 416-3063



Factory Tour



Website